

BY YOUR SIDE **ANYWHERE**



HOME CARE

Home Care aids for disabled and elderly patients with mobility difficulties:

lift chairs, commode chairs, walkers, rollators, crutches, three legs cane, quad cane, canes, supports, IV poles, tables, toilet raisers, shower and bath seats, rehabilitation, devices for the prevention of decubitus, etc.

DUBLINO

Modern design **lift chair with one and two motors**. It is the ideal chair to be adapted to any type of furniture, even to the smallest spaces, to have **maximum relax and stand up without any effort**.



GREY FABRIC
Ref.401000 1 motor
Ref.401010 2 motors



LIGHT BLUE FABRIC
Ref.401001 1 motor
Ref.401011 2 motors



BORDEAUX FABRIC
Ref.401002 1 motor
Ref.401012 2 motors

Class I Medical Device

FEATURES

- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with greater padding for neck and head support
- > Rounded armrests with wide support base
- > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The 1-motor version reclines the backrest and lifts the footrest in a single movement
- > The 2-motors version allows you to separately adjust the inclination of the backrest and the movement of the footrest
- > Maximum backrest inclination 135 °
- > N. 2 rear wheels
- > Handy side pocket for magazines and remote control
- > Gran Foulard in the same fabric and color as the armchair to protect the upholstery

TECHNICAL FEATURES

1-motor version - Electric adjustment in 3 positions



2-motors version - Electric adjustment in 5 positions



Dimensions in cm									
DUBLINO 1 motor	73	106	91	166	49	47	54	35	130
DUBLINO 2 motors	73	106	91	166	49	47	54	40	130

MADRID

Madrid is the **lift chair** with one and two motors which, thanks to its **slim** line, takes up space in width. However, its design does not compromise the comfort guaranteed by the **soft padding** of the backrest and the standard depth of the seat, completely in line with that of other armchairs. Combining **classic and modern** style, Madrid is available in four colours to meet all furnishing needs and tastes.



BORDEAUX MICROFIBER
Ref.401050 1 motor
Ref.401060 2 motors



CORAL MICROFIBER
Ref. 401051 1 motor
Ref. 401061 2 motors



BLUE MICROFIBER
Ref. 401052 1 motor
Ref. 401062 2 motors



OLIVE GREEN MICROFIBER
Ref. 401053 1 motore
Ref. 401063 2 motori

Class I Medical Device

FEATURES

- > Steel and wood frame
 - > Seat padding in non-deformable Polyurethane foam
 - > Ergonomic backrest with extra padding to ensure maximum comfort for people of all stature
 - > Seagull-wing armrests with a wide support base
 - > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The 1-motor version reclines the backrest and lifts the footrest in a single movement
 - > The 2-motor version allows you to separately adjust the inclination of the backrest and the movement of the footrest
 - > Maximum backrest inclination 135 °
 - > N. 2 rear wheels
 - > Handy side pocket for magazines and remote control
 - > Gran Foulard in the same fabric and color as the armchair to protect the upholstery

TECHNICAL FEATURES

1-motor version - Electric adjustment in 3 positions



2-motors version - Electric adjustment in 5 positions



Dimensions in cm									
MADRID 1 motor	72	100	90	165	50	47	50	35	130
MADRID 2 motors	72	100	90	165	50	47	50	40	130

VIENNA

Lift chair with two motors, it has a soft and welcoming high backrest and rest legs joined to the seat, thus ensuring maximum comfort. This electric lift chair is intended both for people with disabilities and for those who simply want to enjoy maximum comfort and relax at home. It has a **Roller System**, wheels system to facilitate the movement of the chair with person on it.



SAND FABRIC
Ref. 401080



ROYAL BLUE FABRIC
Ref. 401081



COFFEE FABRIC
Ref. 401082


Class I Medical Device

FEATURES

- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with extra padding to ensure the maximum comfort for people of all sizes
- > Rounded armrests with large support base
- > "Stand up" function which allows lifting of the person. With this device the person from the sitting position is accompanied up to an almost erect position, facilitating standing
- > The chair is only available in the version with 2 motors: one allows you to adjust the footrest and the other the inclination of the backrest up to 170° reaching the lying position
- > Roller System mechanism for easy movement
- > Handy side pocket for magazine and remote control
- > Gran Foulard in the same fabric and color as the chair to protect the upholstery

TECHNICAL FEATURES

2-motors - Electric adjustment in 6 positions



Dimensions in cm									
VIENNA 2 motors	82	106	92	165	56	49	52	53	170

GINEVRA

Ginevra **lift chair with one and two motors**, with a classic and elegant line that makes it particularly suitable in any type of environment and pleasant for people of all ages. **The extra padding of the backrest** and the **headrest flaps** ensure maximum comfort in all positions. The lift chair has a **Roller System**, a wheel system that allows to move the chair without any effort even with the person seated on it.



COFFEE MICROFIBER
Ref.401020 1 motor
Ref.401030 2 motors



BEIGE MICROFIBER
Ref.401021 1 motor
Ref.401031 2 motors



BLUE MICROFIBER
Ref.401022 1 motor
Ref.401032 2 motors

Class I Medical Device

FEATURES


- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with extra padding to ensure maximum comfort for people of all stature
- > Rounded armrests with a wide support base
- > Headrest on the sides of the backrest for a comfortable support of the head
- > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The 1-motor version lifts the footrest and reclines the backrest up to a maximum of 135 °
- > The 2-motors version allows you to separately adjust the movement of the footrest and the inclination of the backrest up to 170 °, reaching the lying position
- > N. 2 rear wheels (in the version 1 motor only)
- > Handy side pocket for magazines and remote control
- > Gran Foulard in the same fabric and color as the armchair to protect the upholstery
- > Equipped with Roller system mechanism that allows, by pushing with the foot the bar placed behind the chair, to easily move the chair. Once you have entered the mechanism, the lift function is blocked so as to avoid accidental activations

TECHNICAL FEATURES

1-motor version - Electric adjustment in 3 positions



2-motors version - Electric adjustment in 6 positions



Dimensions in cm									
GINEVRA 1 motor	80	107	90	173	54	53	50	46	170
GINEVRA 2 motors	80	107	90	182	54	53	50	51	170

OSLO

Oslo **lift chair** with soft and welcoming outlines is characterized by the presence of many accessories that enrich the design and ease of use. **The movement mechanism** on four independent wheels and the **push handle** allow you to move the chair easily, even with the person seated on it. The **removable headrest flaps** and the **removable armrests** facilitate translation movements from the chair to the bed and vice versa. Equipped with 2 motors; one allows the classic movements "lying down" and "standing up", while the second is intended for the lifting system that allows the user to lift the level of the chair to facilitate the transition from the chair itself to the bed or the wheelchair, eliminating differences in height.



BLACK PU
Ref.401040



IVORY PU
Ref.401041



Class I Medical Device

FEATURES

- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with soft anatomical padding to offer comfort and a correct position
- > Easily removable armrests to facilitate movement and allow the chair to be brought closer to the table
- > Removable headrest flaps for comfortable head support
- > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The chair is only available in the 2-motors version: one allows both the backrest to be reclined and the footrest lifted in a single movement, and the "stand up" position; the other allows you to horizontally lift the base of the chair up to a maximum of 18 cm from the floor
- > Maximum backrest inclination 135 °
- > N. 4 swivel wheels of 10 cm
- > Push handle to easily move the person inside the home
- > Handy side pocket for magazines and remote control

TECHNICAL FEATURES

Electric adjustment in 4 positions



Dimensions in cm									
OSLO 2 motors	74	109	87	173	56	53	54	50	170

OSLO VM

The uniqueness of the characteristics of the Oslo lift chair is enriched, in the VM version, by the **vibromassage and heating functions**. The complex integrated vibromassage mechanism performs a gentle massage that **stimulates blood circulation, helps revitalize muscle tone** e gives a pleasant feeling of **well-being and relax**. The quality of the materials, the different programs with which to massage different parts of the body, as well as the possibility of raising and lowering the seat, backrest and footrest, also taking advantage of the area it heats for **soothe tensions or stiffenings**, make the Oslo VM lift chair a choice of election for those who want to combine relax to self-care. A comfortable and practical piece of furniture not only for the elderly or people with mobility difficulties but also for those users who like to relax at the end of the day, combining **comfort** to the **benefits of a massage**.



Class I Medical Device

FEATURES

- > The Oslo VM relaxation massage lift chair, in addition to Oslo features, combines the two aspects that are most popular and that many want to find: a better quality of life and one pleasant sensation of comfort and well-being throughout the whole body
- > Taking advantage of the massaging swing system, through eight vibrating motors, the chair takes care of stimulating several body parts (back, lumbar, buttocks and legs) with a variable power, ideal for all needs
- > At shoulder vibrations facilitate melting of tension in the cervical area allowing relaxation of the neck and back muscles
- > In the lumbar region, thanks to controlled and continuous use of the vibromassage treatment, you can have significant benefits in the treatment of muscle stress
- > In the buttocks the massage promotes blood circulation e decongests the muscle groups; it is particularly suitable for people who lead a sedentary life
- > In the area of the legs, at the height of the calves, the vibromassage promotes blood circulation, and, if applied while the chair is in a horizontal position, facilitates blood flow in the peripheral areas of the body
- > An easy to use multifunction remote control with LCD screen to program the timer, choose from 10 massage programs and 10 different intensity levels or one manual selection that allows the user to customize the massage optimizing the result
- > Thanks to a heating system positioned in the lumbar area, functional for those who, in addition to massage, also appreciate the soothing action of heat, the Oslo VM lift chair relaxes and instills relaxation and pleasure
- > The Polyurethane cover, easily manageable in terms of daily care, gives practicality and elegance to the lift chair

TECHNICAL FEATURES

Electric adjustment in 4 positions									
Dimensions in cm									
OSLO VM 2 motors	74	109	87	173	56	53	54	50	170

ZERO GRAVITY

Lift chair with a classic and elegant taste that make it suitable not only for people with walking difficulties but also for those looking for a relaxing armchair. This chair has **3 motors** that allow independent adjustment of the footrest and inclination of the backrest, as well as to assume the **Trendelenburg position** or **anti-shock position** to get better perfusion of vital organs such as brain, heart and kidneys.



MELANGE BROWN FABRIC
Ref. 401090



MELANGE BLUE FABRIC
Ref. 401091



Class I Medical Device

FEATURES


- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with extra padding to ensure greater support
- > Rounded armrests with large support base
- > "Stand up" function which allows lifting of the person from the sitting position to an almost erect position, facilitating standing
- > Equipped with 3 independent motors
- > N. 4 wheels (the 2 rear ones with brake)
- > Handy side pocket for magazine and remote control
- > Gran Foulard in the same fabric and color as the chair to protect the upholstery

TECHNICAL FEATURES

3-motors - Electric adjustment in 8 positions



Dimensions in cm									
ZERO GRAVITY 3 motors	82	107	92	180	56	49	50	55	170



KIT AIRPUMP WITH PRESSURE ADJUSTEMENT AND BUBBLE MATTRESS MC 2500M

Ref. 400751

FEATURES

- > Simple installation, use and maintenance
- > Pressure regulation depending to the patient weight
- > Alternating cycle operation to prevent bedsores
- > The airpump is fitted with hooks to be positioned at the foot of the bed
- > Mattress consisting of 130 atoxic PVC cells
- > Equipped with end flaps for a better fixing
- > All parts of the mattress are washable
- > Repair kit included

MATTRESS TECHNICAL FEATURES

- > Dimensions (including end flaps) 280 cm (L) x 90 cm (W) x 6.5 cm (H)
- > Beige colour
- > Weight 2 kg
- > Max weight 130 kg

AIRPUMP TECHNICAL FEATURES

- > Power supply 220V / 50 Hz
- > Air output > 4.5 liters per minute
- > Pressure range 45 - 110 mmHg
- > Cycle time 12 minutes
- > Power 7 W
- > 1A Fuse
- > Dimensions 24.5cm (L) x 12cm (W) x 9.3cm (H)
- > Weight 930 g



LIGHT
AND
SIMPLE
TO USE



HANGING
HOOKS



END
FLAPS



130 CELLS



MAX
WEIGHT
130 KG

SPARE PARTS

Bubble mattress 2500M

Ref. 409780



KIT AIRPUMP WITH PRESSURE ADJUSTEMENT AND INTERCHANGEABLE ELEMENTS MATTRESS MC 5000M

Ref. 400752

FEATURES

- > Simple installation, use and maintenance
- > Pressure regulation depending on the patient's weight
- > Alternating cycle operation to prevent bedsores
- > The airpump is equipped with hooks to be positioned at the foot of the bed
- > 17 interchangeable elements in atoxic Nylon / PVC; the tubes can be positioned individually by means of special buttons
- > The first 3 elements are static and they remain inflated to offer support to the cervical area ensuring maximum comfort to the patient
- > Equipped with CPR (Cardio Pulmonary Resuscitation) valve for rapid deflation of the mattress in case of urgent interventions
- > Waterproof and breathable material which does not cause friction and holds all the elements together for better compactness
- > Equipped with end flaps for better fixing
- > All parts of the mattress are washable

MATTRESS TECHNICAL FEATURES

- > Dimensions 190 cm (L) x 88 cm (W) x 11 cm (H)
- > Blue colour
- > Weight 4.7 kg
- > Max weight 140 kg

AIRPUMP TECHNICAL FEATURES

- > Power supply 220V / 50 Hz
- > Air output 5/6 liters per minute
- > Pressure range 40 - 100 mmHg
- > Cycle time 12 minutes
- > Power 7 W
- > 1A Fuse
- > Dimensions 22.4 cm (L) x 10.8 cm (W) x 10.3 cm (H)
- > Weight 1.4 kg



LIGHT
AND
SIMPLE
TO USE



HANGING
HOOKS



END
FLAPS



WATERPROOF
AND
BREATHABLE
MATERIAL




17
INTERCHANGEABLE
ELEMENTS



PILLOW
FUNCTION



MAX
WEIGHT
140 KG



Ref. 400755

3 SECTIONS STATIC PU FOAM MATTRESS

FEATURES

- > The mattress is divided into three sections to fit one or two functions beds giving to the patients the best comfort in all positions
- > The truncated pyramid design ensures air circulation between the mattress and the patient maintaining the ideal environment
- > Density 30 kg/m³

TECHNICAL FEATURES

- > Dimensions: 190 (3 x 63 cm) x 85 x 12,5 cm
- > Max weight 130 kg

ACCESSORIES

- Polyester and Cotton Fabric Cover with zipper 198 x 90 x 15 cm
- Waterproof PU Cover with zipper 198 x 90 x 15 cm

Ref. 409790

Ref. 409791

MULTI-ZONED STATIC PU FOAM MATTRESS

FEATURES


- > The single section Polyurethane foam mattress is made up of three layers of different densities to ensure differentiated pressure discharge at specific body areas
- > Lateral containment sides density 32 kg/m³
- > Head zone density 29 kg/m³
- > Trunk-pelvis zone density 34 kg/m³
- > Feet perforated zone density 31 kg/m³
- > Equipped with containment cover

TECHNICAL FEATURES

- > Dimensions: 198 x 90 x 15 cm
- > Max weight 130 kg

ACCESSORIES

- Polyester and Cotton Fabric Cover with zipper 198 x 90 x 15 cm
- Waterproof PU Cover with zipper 198 x 90 x 15 cm



Ref. 400756

Ref. 409790

Ref. 409791



The Air cells cushions Medicare have been designed to protect the patient (wheelchair users, elderly people and people who have to stay in the sitting position for a long period) from the risk of forming pressure sores.




AIR CELLS CUSHION
in PVC and Medical rubber

FEATURES

- > The air cells cushion protects against the risk of pressure sores and tissue injuries
- > The combination of materials, PVC and Medical rubber, allows to maintain the cells soft granting comfort and stability to the patient
- > The micro-interchange cells ensure even weight distribution
- > The air cells cushion is highly conformable to the patient's anatomy whose ischemic points are completely immersed in the cushion favoring regular pressure distribution
- > Cushion cover in breathable elastic fabric that allows air to circulate with every movement of the patient
- > The base of the cushion cover in non-slip material secures the cushion to the chair and helps to increase the stability of the user, guaranteeing thus its therapeutic efficacy
- > Easy and adjustable inflation
- > Lightweight and easy to carry and clean
- > Repair kit and inflation pump included



Pressure level before
the use of the cushion



Pressure level after
the use of the cushion

TECHNICAL FEATURES

REF	HEIGHT	VALVS	DIMENSIONS	COLOUR
400760	6 cm	1	40X40 cm	Grey
400761	6 cm	1	45X40 cm	Grey



LIGHT
AND
SIMPLE
TO USE



MAX
WEIGHT
130 KG

VISCOELASTIC ANATOMIC CUSHION


FEATURES

- > Anatomically shaped cushion in compact viscoelastic polymer which guarantees comfort and stability to the patient and prevents the genesis of pressure sores
- > The elastic PU cover, washable and removable by zipper, is fire retardant and impermeable to liquids
- > Non-slip base cover
- > Equipped with velcro strips for fixing the cushion to the wheelchair and avoiding it slipping

TECHNICAL FEATURES

REF	HEIGHT	DIMENSIONS
400785	6 cm	41X41 cm
400786	6 cm	46X41 cm





Ref. 400775

Ref. 400776

CONTOUR CUSHION WITH SELF-MODELING GEL

FEATURES

- > Composite anti-decubitus cushion with preformed anatomical base integrated with self-modeling gel, indicated for patients who use wheelchairs and long-term patients
- > Made up of viscoelastic "slow memory" foam and fluid gel to reduce the pressure on contact areas most at risk of decubitus
- > Removable PU cover, stretch, washable, with base in non-slip material
- > Equipped with handle to facilitate transport

TECHNICAL FEATURES

- > Dimensions 41x41x7 cm
- > Dimensions 44x44x7 cm

CONTOUR CUSHION WITH SELF-MODELING FLUIDS


FEATURES

- > Composite anti-decubitus cushion with preformed anatomical base integrated with self-modeling fluids, suitable to prevent the formation of pressure sores in patients who use wheelchairs and long-term patients who need a comfortable and effective support
- > Made up of a pre-formed base in compact, high density and high resilience foam, and a layer of viscoelastic fluid gel with a cooling effect
- > Removable PU cover, stretch, washable, with base in non-slip material
- > Equipped with handle to facilitate transport

TECHNICAL FEATURES

Ref. 400780 > Dimensions 41x41x7 cm

Ref. 400781 > Dimensions 44x44x7 cm



Ref. 400770

Ref. 400771

FLUID GEL CUSHION

FEATURES

- > The fluid gel cushion guarantees the optimization of pressures to limit the risks of prolonged capillary occlusion in people with reduced mobility
- > The interior is made up of fluid gel bag enclosed between 2 layers of expanded Polyurethane
- > Equipped with stretch, washable, non-slip removable PU cover

TECHNICAL FEATURES

- > Dimensions 40x40x5 cm
- > Dimensions 45x40x5 cm



Ref. 400765

Ref. 400766

SILICONE COTTON WOOL CUSHION

FEATURES

- > Three-section cushion for the prevention of pressure ulcers
- > The silicone cotton wool padding guarantees softness and breathability to the cushion
- > Covered on both sides with 100% sanforized cotton fabric
- > Equipped with double velcro to secure the cushion to the chair and/or to the wheelchair

TECHNICAL FEATURES

- > Dimensions 40x40x10 cm
- > Dimensions 45x40x10 cm

SILICONE COTTON WOOL CUSHION WITH HOLE

FEATURES

- > Cushion with central hole for the prevention of pressure ulcers in particularly sensitive areas and for prostatic problems
- > The silicone cotton wool padding guarantees softness and breathability to the cushion
- > The 100% sanforized cotton fabric cover allows the skin to stay dry and cool

TECHNICAL FEATURES

- Ref. 400769 > Dimensions 42x42x10 cm



Ref. 400791

SILICONE COTTON WOOL HEEL PROTECTOR

FEATURES

- > The heel protection guarantees to the patient the optimization of the pressures on riskful areas of the body
- > The silicone cotton wool padding reduces the heel pressure and avoids friction and consequent skin irritation
- > The 100% sanforized cotton fabric cover allows the skin to be dry and fresh
- > Equipped with a practical adjustable velcro closure

TECHNICAL FEATURES

- > Dimensions 14x10x3 cm

SILICONE COTTON WOOL ELBOW PROTECTOR

FEATURES

- > The elbow protection guarantees to the patient the optimization of pressures on riskful areas of the body
- > The silicone cotton wool padding reduces the elbow pressure and avoids friction and consequent skin irritation
- > The 100% sanforized cotton fabric cover allows the skin to be dry and fresh
- > Equipped with a practical adjustable velcro closure

TECHNICAL FEATURES

- Ref. 400790 > Dimensions 30x15x11 cm



WOODEN CANE

FEATURES

- > Height 93 cm
- > Non-slip rubber tip
- > Available in brown and black



Ref. 400251	Curved handle wooden cane	ø 22mm	brown
Ref. 400252	Curved handle wooden cane	ø 22mm	black
Ref. 400271	Hook handle wooden cane	ø 25mm	brown
Ref. 400272	Hook handle wooden cane	ø 25mm	black
Ref. 400261	T-handle wooden cane	ø 22mm	brown
Ref. 400262	T-handle wooden cane	ø 22mm	black



FOLDABLE ALUMINIUM CANE

FEATURES

- > Anodized aluminium structure
- > 10 levels of height adjustment
- > Ergonomic wooden handle
- > Non-slip tip

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
70-93 cm	0,30 kg	100 kg

Ref. 400235



ADJUSTABLE ALUMINIUM CANE

FEATURES

- > Anodized aluminium structure
- > 10 levels of height adjustment
- > Ergonomic wooden handle
- > Non-slip tip
- > Available in 4 patterns

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
70-93 cm	0,30 kg	100 kg



Ref. 400240	Bon ton
Ref. 400241	Optical
Ref. 400242	Flowers
Ref. 400243	Fantasy



CRUTCHES

FEATURES

- > Structure in light anodized aluminium
- > 10 levels of height adjustment
- > Forearm support
- > Ergonomic handgrips with reflectors on the front
- > Non-slip rubber tip
- > Available in different colours

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
90-112 cm	0,40 kg	120 kg

Ref. 400201	grey
Ref. 400202	blue
Ref. 400203	red
Ref. 400204	green



CRUTCHES with soft handle

FEATURES

- > Structure in light anodized aluminium
- > 10 levels of height adjustment every 2,5 cm
- > Forearm support
- > Ergonomic, non-slip and soft handgrips with reflectors on the front
- > Anti-noise ring on the sliding tube
- > Non-slip rubber tip
- > Available in grey and black colours

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
92-116 cm	0,45 kg	120 kg

Ref. 400211	gray
Ref. 400212	black





Ref. 400230 **boy**
Ref. 400231 **adult**

UNDERARM CRUTCHES

FEATURES

- > Structure in light anodized aluminium
- > 9 levels of height adjustment
- > Ergonomic forearm support and handle
- > Non-slip rubber tip

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
97-114 cm	0,80kg	110 kg
113-134 cm	0,90kg	110 kg



THREE LEGS CANE

FEATURES

- > Structure in chromed aluminium
- > Height adjustable
- > Ergonomic handle
- > Non-slip rubber tip
- > Reinforced base

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
74-97 cm	0,76 kg	113 kg

Ref. 400301

QUAD CANE

FEATURES

- > Structure in chromed aluminium
- > Height adjustable
- > Ergonomic handle
- > Non-slip rubber tips
- > Reinforced base

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
74-97 cm	0,78 kg	113 kg

Ref. 400351



HD QUAD CANE

FEATURES

- > Wide base structure in anodized aluminium
- > Height adjustable
- > Heavy duty
- > Ergonomic handgrips
- > Non-slip rubber tips

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
68-91 cm	0,8 kg	136 kg

Ref. 400355 **green**
Ref. 400356 **bordeaux**



Ref. 400601

WALKER reciprocating

FEATURES

- > Aluminium reciprocating walker
- > Foldable by one button
- > Ergonomic handgrips
- > 8 levels height adjustment
- > 4 non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
78-96 cm	47 X 46 cm	2,4 kg	136 kg



WALKER foldable fix

FEATURES

- > Aluminium fix walker
- > Foldable by two buttons
- > Ergonomic handgrips
- > 8 levels height adjustment
- > 4 non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
79-96 cm	48 X 44 cm	2,6 kg	136 kg

Ref. 400501

WALKER double handgrips

FEATURES

- > Aluminium fix walker
- > Foldable by one button
- > 2 height handgrips to ease the patient to stand up
- > Ergonomic handgrips
- > Height adjustable
- > 4 non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
74-87 cm	56 X 52 cm	2,7 kg	136 kg



Ref. 400510



WALKER 2 wheels - 2 tips

FEATURES

- > 2 wheels and 2 tips aluminium walker
- > Foldable by one button
- > Ergonomic handgrips
- > 8 levels height adjustment
- > 2 non-slip rubber tips and 2 x 5" wheels

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
79-95 cm	61 X 52 cm	3 kg	136 kg

Ref. 400551

WALKER
2 castors - 2 wheels with brake

FEATURES

- > 2 x 3" castors and 2 x 3" wheels with brake aluminum walker
- > Foldable by one button
- > Ergonomic handgrips
- > Height adjustable

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
83-100 cm	61 X 55 cm	3 kg	110 kg

Ref. 400552



HIGH WALKER
with underarm supports

FEATURES

- > Foldable high walker with height adjustable soft PU underarm supports
- > Steel frame with 4 castors (2 of them with brake)
- > Equipped with padded seat

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
100-140 cm	78 X 55 cm	60 cm	15 kg	120 kg

Ref. 400402



Ref. 400401

HIGH WALKER
with table

FEATURES

- > Foldable high walker with padded table and handles
- > Steel frame with 4 x 100mm castors (2 of them with brake)
- > Height adjustable table and handles
- > Equipped with padded seat, backrest and object holder hook

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
110-135 cm	73 X 64 cm	55 cm	11 kg	136 kg

WALKER
ACCESSORIES

5" CASTOR EXTENSION LEGS

Ref. 409831



3" CASTOR EXTENSION LEGS

Ref. 409830



5" WHEEL EXTENSION LEGS

Ref. 409833



3" WHEEL EXTENSION LEGS

Ref. 409832



PRESSURE BRAKE
EXTENSION LEGS

Ref. 409835



PRESSURE BRAKE AND
3" WHEELS EXTENSION LEGS

Ref. 409834



ROLLATOR foldable 2 wheels

FEATURES

- > Foldable steel frame rollator with 2 x 5" wheels and 2 tips
- > Ergonomic handgrips
- > Height adjustable
- > Padded seat

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
84-94 cm	61 X 55 cm	55-65 cm	6 kg	120 kg

Ref. 400605



Ref. 400606

ROLLATOR SKINNY foldable 2 wheels

FEATURES

- > Foldable aluminium frame rollator
- > 2 x 6" front wheels and 2 back tips
- > Easily foldable for storage and transport
- > Ergonomic handgrips
- > Padded backrest band
- > Height adjustable in 5 different positions
- > Removable shopping bag

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
78-97 cm	60 X 56cm	47-54 cm	5.5 Kg	136 Kg

ROLLATOR foldable 4 wheels with seat

FEATURES

- > Foldable aluminium frame rollator with 8" front castors and back wheels
- > Handle brake on the back wheels
- > Ergonomic handgrips
- > Height adjustable handles
- > Padded seat and backrest band
- > Equipped with basket and cane holder

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
81-92 cm	61 X 61 cm	54 cm	7.5 kg	120 kg

Ref. 400610



ROLLATOR CUBE
foldable and compact

FEATURES

- > Elegant, light, foldable aluminium frame rollator
- > Oval frame completely foldable to grant minimum storage dimensions (63x30x35 cm) and to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Handle brake
- > Height adjustable ergonomic handgrips for extra comfort
- > Equipped with seat, backrest, shopping bag with reflector and cane holder

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
80-93 cm	63X68 cm	54 cm	7,5 kg	136 kg



Ref. 400620



Ref. 400615

ROLLATOR ECO-CUBE
foldable and compact

FEATURES

- > Elegant foldable rollator which combines the lightness of aluminium frame to the refinement of the wood-effect coloring
- > Oval frame completely foldable to grant minimum storage dimensions (90x26,5x33 cm) and to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Handle brake
- > Height adjustable ergonomic handgrips for extra comfort
- > Equipped with seat, backrest, and cane holder
- > Shopping bag with reflector

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
78-93 cm	69X61,5 cm	55,5 cm	5,5 Kg	136 kg





Ref. 400625

ROLLATOR URBAN 2 functions 4 wheels

FEATURES

- > Practical and comfortable aluminium frame rollator convertible into transit wheelchair
- > Double crossbar foldable frame to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Handle brake
- > Height adjustable ergonomic handles
- > Equipped with seat, backrest, shopping bag with reflector and cane holder
- > Flip-up footrests with heel strip

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
84-97 cm	63X75 cm	55.5 cm	9.2 Kg	136 kg

ROLLATOR URBAN AB with PU forearm pads

FEATURES

- > Height adjustable (5 positions) aluminium frame rollator
- > Side-to-side foldable frame to grant minimum storage dimensions (23x61x100 cm) and to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Height adjustable soft ergonomic forearm pads
- > Handle steel brakes
- > Equipped with fabric seat and cane holder
- > Detachable shopping bag with reflector

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
100-115 cm	70X61 cm	60.5 cm	7.6 Kg	136 Kg

Ref. 400626



SPARE PARTS

Ref. 409825	Wheels for rollator Ref. 400551	Ref. 409915	Rubber tips for Ref. 400301
Ref. 409827	100 mm back wheels with brake for rollator Ref. 400402 (pair)	Ref. 409920	Rubber tips for Ref. 400351
Ref. 409828	100 mm front wheels for rollator Ref. 400402 (pair)	Ref. 409925	Rubber tips for Ref. 400201-202-203-204
Ref. 409900	Rubber tips for Ref. 400551	Ref. 409926	Rubber tips for ø 22mm wooden canes
Ref. 409905	Rubber tips for Ref. 400601	Ref. 409927	Rubber tips for ø 25mm wooden canes
Ref. 409910	Rubber tips for Ref. 400501	Ref. 409930	PU handgrip for Ref. 400501-400551-400601

COMMODE CHAIR WITH RECLINING BACKREST self-propelling

FEATURES

- > Steel frame
- > Reclining and removable padded backrest, 56 cm high
- > Padded seat with anatomic opening for toilet
- > Made up of washable and fire retardant material
- > 4 PU wheels, 2 x 24" (60 cm) and 2 x 8" (20 cm)
- > Brake on the rear wheels
- > Padded, flip-up and removable armrests
- > Height adjustable elevating footrests with padded legrests
- > Flip-up and lenght adjustable footrests
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
66 cm	43 cm	50 cm	106 cm	25 kg	130 kg

Ref. 400050



HEADREST for commode chairs



Ref. 409730



Ref. 400053

COMMODE CHAIR WITH RECLINING BACKREST 4" wheels

FEATURES

- > Steel frame
- > Reclining and removable padded backrest, 56 cm high
- > Padded seat with anatomic opening for toilet
- > Made up of washable and fire retardant material
- > 4 x 4" swivel wheels, the 2 rear ones with brake
- > Padded, flip-up and removable armrests
- > Height adjustable elevating footrests with padded legrests
- > Flip-up and lenght adjustable footrests
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
56 cm	43 cm	50 cm	106 cm	21 Kg	130 Kg

COMMODE CHAIR WITH RECLINING BACKREST 8" wheels

FEATURES

- > Steel frame
- > Reclining and removable padded backrest, 56 cm high
- > Padded seat with anatomic opening for toilet
- > Made up of washable and fire retardant material
- > 4 x 8" swivel wheels, the 2 front ones with brake
- > Padded, flip-up and removable armrests
- > Height adjustable elevating footrests with padded legrests
- > Flip-up and length adjustable footrests
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
62 cm	43 cm	50 cm	106 cm	22 kg	130 Kg



Ref. 400055



COMMODE CHAIR

FEATURES

- > Height adjustable steel frame
- > Flip-up armrests
- > Padded seat
- > Water-repellent leather seat
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
55 cm	43.5 cm	50-57 cm	95-103 cm	10 kg	136 kg

Ref. 400034

COMMODE CHAIR with wheels

FEATURES

- > Steel frame with 4 x 5" swivel wheels, 2 of which with brake
- > Flip-up armrests
- > Adjustable and removable footrests
- > Padded seat
- > Water-repellent leather seat
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
60 cm	43 cm	50 cm	88 cm	12 kg	120 kg



Ref. 400035



3 IN 1 COMMODE CHAIR

FEATURES

- > 3 uses in 1: toilet raiser, commode chair, toilet seat
- > Steel frame
- > Height adjustable
- > Plastic armrest and round bucket with lid
- > Non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
57-72 cm	47 X 58 cm	37-52 cm	4,8 Kg	136 kg

Ref. 400030

SHOWER AND COMMODE CHAIR self-propelling

FEATURES

- > Reinforced anti-rust and anti-corrosion aluminium frame
- > Lightweight and versatile: it can be used as a toilet seat
- > PU seat with anatomic opening for toilet to allow the use of the bucket
- > 2 x 24" PU wheels and 2 x 5" wheels with brake
- > Removable PU backrest
- > Flip-up armrests
- > Swiveling, removable footrests with height-adjustable footplates

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
64 cm	46 cm	50 cm	99 cm	18 kg	120 kg

Ref. 400040



SHOWER AND COMMODE CHAIR 5" wheels

FEATURES

- > Reinforced anti-rust and anti-corrosion aluminum frame
- > Lightweight and versatile: it can be used as a toilet seat
- > PU seat with anatomic opening for toilet to allow the use of the bucket
- > 4 x 5" PU wheels with brake
- > Removable PU backrest
- > Flip-up armrests
- > Swiveling, removable footrests with height-adjustable footplates

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
55 cm	46 cm	50 cm	99 cm	15 kg	120 kg

Ref. 400045





Ref. 400101

SHOWER CHAIR

FEATURES

- > Aluminium frame
- > Ergonomic seat with holes for water drainage
- > Side handles
- > Height adjustment at 8 different levels
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	SEAT HEIGHT	WEIGHT	MAX WEIGHT
30X50 cm	37-50 cm	1.8 kg	113 Kg

WALL SHOWER SEAT

FEATURES

- > Flip-up aluminium frame to be fixed to the wall
- > Floor legs
- > Ergonomic seat with holes for water drainage
- > Side handles
- > Height adjustable
- > Materials resistant to corrosion

SEAT	SEAT HEIGHT	WEIGHT	MAX WEIGHT
30X50 cm	40-55 cm	4.8 kg	113 Kg



Ref. 400102



Ref. 400110

SHOWER CHAIR WITH BACKREST

FEATURES

- > Aluminium frame
- > Ergonomic seat and backrest with holes for water drainage
- > Side handles
- > Height adjustment at 8 different levels
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	BACK	SEAT HEIGHT	BACK HEIGHT	WEIGHT	MAX WEIGHT
30X50 cm	35X50 cm	37-50 cm	70-87 cm	2.7 kg	120 kg

FOLDABLE SHOWER CHAIR

FEATURES

- > Aluminum foldable frame
- > Ergonomic backrest and seat with holes for water drainage
- > Height adjustable
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	BACK	SEAT HEIGHT	BACK HEIGHT	WEIGHT	MAX WEIGHT
40x41 cm	20x48 cm	46-54 cm	88-96 cm	8,2 Kg	113 Kg

Ref. 400109



Ref. 400115

SWIVEL SHOWER CHAIR

FEATURES

- > Steel frame
- > Ergonomic swivel PVC seat
- > Armrests for more comfort

SEAT	BACK	BASE	WEIGHT	MAX WEIGHT
40X37 cm	40X37 cm	73X52 cm	4,3 kg	120 kg

ROUND SHOWER SEAT

FEATURES

- > Aluminium frame
- > Round seat
- > Height adjustment at 8 different levels
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	SEAT HEIGHT	WEIGHT	MAX WEIGHT
Ø32 cm	38-55 cm	1,6 kg	100 kg

Ref. 400100



BATH BOARD WITH HANDLE AND ROTATING SEAT

FEATURES

- > High resistance PVC seat
- > Non-slip rotating seat
- > Safety handle and adjustable brackets for a better fitting to the bathtub's dimensions

BASE	WEIGHT	MAX WEIGHT
75X36 cm	1,5 kg	120 kg



Ref. 400125



Ref. 400130

SAFETY HANDLE FOR BATHTUB

FEATURES

- > Steel safety handle to facilitate entry and exit from the bathtub
- > Fixing clamp to the bathtub (maximum width 17.8 cm)

WEIGHT	MAX WEIGHT
2,8 kg	120 kg

TRANSFER SEAT FOR BATHTUB

FEATURES

- > Aluminium frame
- > Non-slip PVC seat
- > Backrest with handle
- > Height adjustable
- > Non-slip suction cup tips
- > Materials resistant to corrosion

SEAT	BACK	SEAT HEIGHT	BACK HEIGHT	WEIGHT	MAX WEIGHT
70X41 cm	35X45 cm	46-56 cm	81-91 cm	4,4 kg	120 kg



Ref. 400120

TOILET RAISER

FEATURES

- > Ergonomic seat for more comfort
- > Washable PVC material with high resistance
- > Front and rear opening to ease accessibility for personal hygiene
- > Fixing system with adjustable side brackets



Ref. 400020 - H 10 cm
Ref. 400025 - H 15 cm



Ref. 400010 - H 10 cm
Ref. 400015 - H 15 cm

TOILET RAISER with lid

FEATURES

- > Ergonomic seat for more comfort
- > Washable PVC material with high resistance, with lid
- > Front and rear opening to ease accessibility for personal hygiene
- > Fixing system with adjustable side brackets



Ref. 400005

TOILET RAISER with armrests

FEATURES

- > Ergonomic seat
- > Simple and safe front fixing system
- > It raises the toilet level of 5"
- > Washable PVC material
- > Front opening to ease accessibility for personal hygiene



Ref. 400028

SOFT TOILET RAISER H 10 cm

FEATURES

- > Ergonomic seat for more comfort
- > Washable soft PU material with high resistance
- > Front opening to ease accessibility for personal hygiene
- > Fixing system with adjustable side brackets



IV STAND

FEATURES

- > Steel frame
- > Height adjustable
- > Equipped with 4 hooks
- > Base with 5 castors

HEIGHT ADJUSTABLE	MAX WEIGHT
133-223 cm	2kgxhook

Ref. 400651



Ref. 400701

BACKREST

FEATURES

- > Bed backrest
- > Steel frame
- > Adjustable in 5 different positions
- > Made of breathable nylon with padded headrest
- > Blue color

BASE	HEIGHT ADJUSTABLE	MAX WEIGHT
54x58 cm	30-54 cm	120 kg

EXAMINATION COUCH

FEATURES

- > Steel frame 190(L) x 65(W) x 75(H) cm
- > PU foam padding and washable leather
- > Adjustable in 5 different positions
- > Roll holder
- > Adjustable backrest

MAX WEIGHT
180 kg

Ref. 400660





Ref. 400951

OVERBED TABLE double tilting desk top

FEATURES

- > Steel frame
- > Double tilting desk top surrounded by edges
- > Height adjustable
- > With 4 swivel wheels, 2 of which with brake
- > Dark brown color

SMALL DESK TOP	BIG DESK TOP	HEIGHT ADJUSTABLE
40 X 20 cm	40 X 68 cm	73-125 cm

OVERBED TABLE tilting desk top

FEATURES

- > Steel frame
- > Tilting desk top
- > Height adjustable
- > With 4 swivel wheels, 2 of which with brake
- > Brown color

DESK TOP HEIGHT ADJUSTABLE	HEIGHT ADJUSTABLE
41 X 65 cm	71-110 cm

Ref. 400952



Ref. 400955

TILTING OVERBED TABLE pneumatic adjustment

FEATURES

- > Steel frame
- > Tilting (0° to 330°) desk top, washable, fire retardant, scratch resistant and highly durable
- > Pneumatic height adjustment
- > With 4 swivel wheels, 2 of which with brake
- > Gray color

TABLE TOP HEIGHT ADJUSTABLE	HEIGHT ADJUSTABLE
48 X 78 cm	76-117 cm

PEDAL EXERCISER

FEATURES

- > Steel frame
- > Strength adjustable pedal
- > Triangle anti-slips tips

Ref.400630



Ref.400631

PEDAL EXERCISER digital display

FEATURES

- > Steel frame
- > Digital display for timer, counting number of repetitions, speed and consumed calories
- > Triangle non-slip tips for more stability and safety

SWIVEL STOOL

FEATURES

- > Plastic base
- > Padded PU leather seat
- > Height adjustable

Ref.400668



Ref.400665

STEPPER

FEATURES

- > Steel frame
- > 2 steps
- > Non-slip footrests
- > Non-slip tips