

WHEN A VALID SUPPORT IS NEEDED



SUPPORTS AND BRACES

The partnership with **Conwell Medical Co. Ltd.**, a company specialized in design and manufacturing of a wide range of support and braces, is the result of a keen research. The synergy between the two companies results into a **complete range of patented and professional solutions** for the body parts at risk of injury.

CONWELL Medical is a fast growing company specialized in the field of orthopedic braces and supports. Founded in Taiwan in 2005, Conwell Medical Co., Ltd. is a manufacturer, distributor and exporter of its patented brand Conwell. Conwell offers a complete range of orthopedic products for the major injury-prone areas of the body.


"Be Well with Conwell" is the mission of the Company. To accomplish that, Conwell rigorously selects the materials and constantly innovate new designs to ensure that the products provide effectiveness and comfort to customers.

Product quality control is a priority.

In 2009, Conwell received the quality standard ISO 9001 and ISO 13485.

**To keep its products to a level
of excellence CONWELL has chosen
COOLMAX® and COOLPLUS®**

MOISTURE MANAGEMENT + BREATHABILITY = TOTAL COMFORT



COOLMAX® fabric


Why choosing COOLMAX® fabric?

For the simple reason of staying dry and comfortable when you start to perspire. Sweating is one of the body's natural mechanisms to stay cool. Engineered with special technology to wick moisture and enhance drying, COOLMAX® fabric helps you stay dry and energized - exactly what you need to feel comfortable not only during the day, but also through your toughest workouts.

Today COOLMAX® fabric appeals to the fashion conscious who prefer a cool, soft and high performance fabric that fits every lifestyle. COOLMAX® fabric is the confident, sophisticated choice for active consumers.

- COOLMAX® fabric pulls moisture away from your skin.
- COOLMAX® fabric absorbs and spreads moisture out across fabric to enhance Moisture Vapor Transmission rate.
- Air moves in to keep your body cool and dry.

COOLMAX® is a trademark of INVISTA.



COOLPLUS® fabric

Why choosing COOLPLUS® fabric?

It keeps you cool and dry.

COOLPLUS® is a new fiber that features water diffusion and wicking properties. These fibers have excellent texture yielded several of micro-slits within single fiber. As a result the fabric does not stick to the skin and shows cool crispy wearing comfortable due to its capillary action.

In addition, the dry touch can be maintained even when the part is still wet.

INDICATIONS

- > Mild sprains of the cervical spine
- > Stiff neck
- > Cervical arthrosis

FEATURES

- > Made of polyurethane foam
- > Covered with stockinette for maximum patient comfort
- > Hook and loop closure for ease of adjustment
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
S	36 cm – 42 cm	7 cm
M	39 cm – 45 cm	7 cm
L	41 cm – 47 cm	9 cm
XL	44 cm – 50 cm	9 cm

**SOFT
CERVICAL
COLLAR****5101****PLASTIC
ADJUSTABLE
COLLAR**
5103**INDICATIONS**

- > Support for minor sprains of the cervical spine
- > Stiff neck
- > Discopathies
- > Cervical arthrosis
- > Post-operative
- > Vertiginous crises of cervical origin

FEATURES

- > Semi-rigid polypropylene structure
- > Foam padding with ventilation holes for comfort of the patient
- > Adjustable in height and circumference
- > Hook and loops closure for easy adjustment
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
S	29 cm – 34 cm	6 cm
M	34 cm – 39 cm	7 cm
L	39 cm – 44 cm	9 cm
XL	44 cm – 49 cm	10 cm

INDICATIONS

- > Support for minor sprains of the cervical spine
- > Stiff neck
- > Discopathies
- > Cervical arthrosis
- > Post-operative
- > Vertiginous crises of cervical origin

FEATURES

- > Semi-rigid polypropylene structure
- > Foam padding with ventilation holes for comfort of the patient
- > Adjustable in height and circumference
- > Hook and loops closure for easy adjustment
- > Special design for excellent support
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
S	29 cm – 34 cm	6 cm
M	34 cm – 39 cm	7 cm
L	39 cm – 44 cm	9 cm

**PLASTIC
ADJUSTABLE
COLLAR
with chin
support****5104**

CERVICAL COLLAR with trachea opening h. 4 1/4"

5108



CERVICAL COLLAR with trachea opening h. 5 1/4"

5109

INDICATIONS

- > Support for minor cervical sprains
- > Stiff neck
- > Discopathies
- > Cervical arthrosis
- > Post-operative
- > Vertiginous crises of cervical origin
- > Useful in case of tracheotomy

FEATURES

- > Lightweight and compact bivalve structure
- > Large trachea opening for easy emergency tracheotomies
- > Hook and loop side closures for an easy fit
- > X-ray transparent
- > Available in different heights

SIZE	NECK CIRCUMFERENCE
S	25 cm – 33 cm
M	33 cm – 40 cm
L	40 cm – 48 cm
XL	48+ cm

DELUXE ADJUSTABLE BIVALVE CERVICAL COLLAR

51100

INDICATIONS

- > Cervical trauma
- > Cervical hernias
- > Instability of the cervical tract due to neurological pathologies
- > Cervicoarthrosis
- > Discoarthrosis
- > Myogenic stiff necks
- > Post-operative immobilization

FEATURES

- > Bivalve design made with light and quality materials to fit any morphological profile and improve patient compliance especially in cases of prolonged use
- > Simple height adjustment: this allows the collar to vary the traction on the patient's cervical area
- > The sternal support of the collar increases the stability of the device
- > Wide trachea opening for easy emergency tracheotomies
- > Soft foam padding prevents pressure sores and grants comfort
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
ONE SIZE	21 cm – 49 cm	Adjustable



CLAVICLE BRACE

5201

INDICATIONS

- > Immobilization of the clavicle
- > Subluxation of the joint acromio-clavicular

FEATURES

- > Soft foam padding for patient comfort
- > Padded closed cell foam axillary straps covered with breathable stockinette
- > Hook and loop closure for ease of adjustment

SIZE	CHEST CIRCUMFERENCE
XS	44 cm - 54 cm
S	54 cm - 64 cm
M	64 cm - 76 cm
L	76 cm - 89 cm
XL	89 cm - 104 cm



SHOULDER SUPPORT

5202



INDICATIONS

- > Immobilization of the shoulder
- > Sprains
- > Tendinitis
- > Bursitis

FEATURES

- > Soft foam padded shoulder straps for patient comfort
- > Length-adjustable
- > Ambidextrous

SIZE	ARM CIRCUMFERENCE
S/M	23 cm – 32 cm
L/XL	31 cm – 40 cm

DELUXE SHOULDER IMMOBILIZER

52070

INDICATIONS

- > Immobilization of the shoulder
- > Dislocations
- > Subluxations
- > Post-surgery immobilization

FEATURES

- > Soft foam padded bands to ensure comfort to the patient
- > Two over-the-shoulder criss-cross straps to increase joints stability
- > Length-adjustable straps
- > Final stabilization band with hook and loop closure
- > Hook and loop closure closures for an easy fit
- > Ambidextrous

SIZE

ONE SIZE



SHOULDER IMMOBILIZER 5206

INDICATIONS

- > Treatment of dislocations
- > Post-surgery immobilization

FEATURES

- > Elastic waste band construction for compressive support
- > Durable wrist and humeral cuff with contact closure for an easy fit
- > Ambidextrous

SIZE	WAIST CIRCUMFERENCE
S/M	66 cm – 89 cm
L/XL	89 cm – 112 cm



MESH ARM SLING 5204

INDICATIONS

- > Immobilization and support for mild spains and strains
- > Dislocations of the shoulder complex
- > Post fractures

FEATURES

- > Envelope design with thumb support
- > Hook and loop closure for ease of adjustment
- > Ambidextrous

SIZE	ARM LENGHT (from elbow to wrist)
S	28 cm – 33 cm
M	33 cm – 38 cm
L	38 cm – 43 cm
XL	43 cm – 48 cm



DELUXE ARM SLING 5205

INDICATIONS

- > Immobilization and support for mild spains and strains of the shoulder complex
- > Bruises and bursitis
- > Immobilization for fractures of the shoulder

FEATURES

- > Envelope design with thumb support
- > Shoulder lenght adjustable strap with foam pad for great comfort
- > Hook and loop closure for ease of adjustment
- > Ambidextrous

SIZE	ARM LENGHT (from elbow to wrist)
S	28 cm – 33 cm
M	33 cm – 38 cm
L	38 cm – 43 cm
XL	43 cm – 48 cm

INDICATIONS

- > Back pain
- > Lumbar / abdominal support

FEATURES

- > High quality breathable elastic fabric construction
- > Soft foam padding offers excellent comfort
- > 7.5 cm reinforced elastic cross-bands that can be tensioned for extra support
- > Hook and loop closure to adjust the compression to the desired intensity

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm



SOFT BACK SUPPORT

5501



BACK SUPPORT

8" 5502

**INDICATIONS**

- > Back pain
- > Lumbosacralgia
- > Muscles contractures
- > Post surgery lumbar support
- > Dorsal-lumbar disc herniation

FEATURES

- > Structure in breathable and resistant mesh elastic fabric
- > 2 removable easily shapable metal back stays integrated into the band to support the lumbar area
- > Tensionable cross elastic bands of reinforcement with hook and loop closure to adjust the compression to the desired intensity
- > Front hook and loop closure

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm

**INDICATIONS**

- > Back pain
- > Lumbar trauma
- > Muscles contractures
- > Post surgery lumbar support

FEATURES

- > Structure in breathable and resistant mesh elastic fabric
- > 4 flexible stays in the back and sides for extra lumbar support
- > Tensionable cross elastic bands of reinforcement with hook and loop closure to adjust the compression to the desired intensity
- > Front hook and loop closure

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm



SACRO LUMBAR SUPPORT

10" 5503

10"

5503



INDICATIONS

- > Back pain
- > Lumbar trauma
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Structure in breathable and resistant mesh elastic fabric
- > 4 flexible stays in the back and sides for extra support
- > Extra EVA foam pad at the back for compression and lumbar stabilization
- > Tensionable cross elastic bands of reinforcement to adjust compression and support
- > Front hook and loop closure

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm



SACRO LUMBAR SUPPORT with back pad

12"

5504



SPINAL BRACE with back pad

5505

INDICATIONS

- > Back pain
- > Lumbar trauma
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Made of breathable and resistant mesh elastic fabric for a better distribution of compression
- > 2 removable easily shapable metal back stays integrated into the band
- > EVA foam pad at the back for lumbar stabilization
- > Adjustable shoulder straps with hook and loop closure to promote good back posture
- > Tensionable cross elastic back bands to adjust lumbar compression
- > Front hook and loop closure for an easy fit

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm



DELUXE SACRO LUMBAR SUPPORT 5506



INDICATIONS

- > Back pain
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Made of breathable high quality elastic fabric
- > Abdominal panels to increase abdominal compression
- > 4 flexible stays at the back for extra support
- > 3" elastic tensionable straps for more stability
- > Hook and loop closure for an easy fit

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm

INDICATIONS

- > Post surgery
- > Post pregnancy
- > Abdominal support

FEATURES

- > Made of high quality mesh elastic fabric
- > High abdominal compression
- > Highly breathable to promote the healing process
- > Front hook and loop closure to ease adjustment

SIZE	WAIST CIRCUMFERENCE
S/M	63 cm – 96 cm
L/XL	84 cm – 117 cm

ABDOMINAL SUPPORT BAND 11"

11"

5517



ABDOMINAL SUPPORT BAND 12"

12"

5508

INDICATIONS

- > Abdominal support
- > Post surgery
- > Post pregnancy

FEATURES

- > High quality durable elastic fabric
- > Made of 4 panels for a conforming fit and adjustable compression
- > Front hook and loop closure to ease adjustment
- > 12" height
- > Useful in case of particularly large surgical incisions

SIZE	WAIST CIRCUMFERENCE
S/M	63 cm – 96 cm
L/XL	84 cm – 117 cm





HYPEREXTENSION SPINE BRACE (SHORT LENGTH)

5514

HYPEREXTENSION SPINE BRACE (REGULAR LENGTH)

5515

INDICATIONS

- > Compression of spinal fractures
- > Vertebral collapse
- > Back lumbar arthrosis
- > Traumatic or pathological symptoms resulting from osteoporosis

FEATURES

- > Three-point light aluminum construction with soft foam padding for more comfort
- > Immobilization of the spine thanks to the three support points on the patient: the sternal plate on the sternum, the dorsal plate on the spine and the pelvic band on the pelvis
- > Adjustable and tilting axillary joints allow correct fitting and maximum functionality
- > Ergonomic safety closure for controlled compression

SIZE	HIPS CIRCUMFERENCE	5514 HEIGHT	5515 HEIGHT
S	60 cm – 75 cm	36 cm – 44 cm	39 cm – 47 cm
M	75 cm – 90 cm	39 cm – 47 cm	42 cm – 50 cm
L	90 cm – 105 cm	42 cm – 50 cm	45 cm – 53 cm
XL	105 cm – 115 cm	45 cm – 53 cm	48 cm – 56 cm



INDICATIONS

- > Back pain
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Made of breathable high quality mesh elastic fabric
- > Abdominal panels to increase abdominal compression
- > 4 flexible stays at the back for extra support
- > 3" elastic tensionable straps for more stability
- > Hook and loop closure for an easy fit

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm



BACK SUPPORT

9"

5519



INDICATIONS

- > Support for elbow instability
- > Pre and post elbow surgery

FEATURES

- > Elbow brace allows controlled range of motion for both extension and flexion
- > Detachable hand rest
- > Detachable shoulder sling
- > Ambidextrous

SIZE	FOREARM CIRCUMFERENCE
S/M	Less than 21 cm
L/XL	22 cm – 30 cm

ROM ELBOW BRACE**5309****ELBOW SUPPORT COOLMAX®****53200****INDICATIONS**

- > Supports the elbow joint area
- > Protects against bruising or scraping

FEATURES

- > Made of COOLMAX® fabric
- > Neoprene can be applied even in case of swollen or injured elbow
- > Hoop and loop closure for an easy fit
- > Ambidextrous

SIZE

ONE SIZE**INDICATIONS**

- > Immobilization and protection for wrist sprains and fractures

FEATURES

- > Anatomic plastic splints above and below the palm to provide excellent immobilization
- > The lightness of the splint together with the high quality soft foam grant an optimal comfort
- > Adjustable wrist wrap-around straps for ease of application
- > Ambidextrous

SIZE

ONE SIZE**WRIST SPLINT****53210 - left****53220 - right**

WRIST SPLINT

9"

5314 - left
5314 - right



INDICATIONS


- > Immobilization and protection of the wrist in case of sprains, tendonitis and arthrosis
- > Post fractures treatment
- > Fractures of the extremity of the radius and ulna
- > Substitute of casting

FEATURES

- > Breathable fabric
- > Adjustable hook and loop closure to offer greater stability
- > Removable aluminum splints to improve immobilization
- > Right and left

SIZE	WRIST CIRCUMFERENCE (left / right)
S	13 cm – 15 cm
M	15 cm – 17 cm
L	17 cm – 19 cm
XL	19 cm – 21 cm

BREATHABLE



INDICATIONS


- > Immobilization and protection of the wrist in case of sprains, tendonitis and arthrosis
- > Post fractures
- > Lesions of the ulnar area of the metacarpophalangeal joint of the thumb

FEATURES

- > Breathable fabric
- > Adjustable hook and loop closures to offer greater stability
- > Removable front and side aluminum splints to improve wrist and thumb immobilization
- > Right and left

SIZE	WRIST CIRCUMFERENCE (left / right)
S	13 cm – 15 cm
M	15 cm – 17 cm
L	17 cm – 19 cm
XL	19 cm – 21 cm

BREATHABLE



THUMB WITH WRIST SPLINT

5315 - left
5315 - right



THUMB WRIST SUPPORT

53030



INDICATIONS

- > Wrist and thumb immobilization and support
- > Joint instability
- > Sprains
- > Arthrosis

FEATURES

- > Made of elastic fabric with thumb support
- > Metal stay insert to keep the thumb in position
- > Adjustable wrist and thumb wrap-around straps to offer extra stability
- > Hook and loop closures for ease of application
- > Ambidextrous

SIZE	WRIST CIRCUMFERENCE
ONE SIZE	14 cm – 28 cm

INDICATIONS

- > Compression of tendon insertions in the elbow
- > Resolution of pain caused by sports injuries and by elbow overexertion (tennis elbow)
- > Prevention of epicondylitis

FEATURES

- > Anatomically shaped neoprene structure
- > Integrated pressure pad
- > Hook and loop closure for compression control
- > Elastic finishing edge

SIZE	ELBOW CIRCUMFERENCE
ONE SIZE	18 cm - 36 cm

TENNIS ELBOW SUPPORT**53070****THUMB WRIST SPLINT****53160****INDICATIONS**

- > Immobilization and support of the thumb
- > Tendinitis
- > Rhizoarthrosis
- > Post fracture

FEATURES

- > Durable and breathable perforated fabric to ensure excellent immobilization and patient comfort
- > Adjustable hook and loop closure to offer greater stability
- > Metal splint to keep thumb in place
- > Ambidextrous

SIZE	WRIST CIRCUMFERENCE
ONE SIZE	14 cm - 28 cm

INDICATIONS


- > Sprains of the thumb
- > Treatment of De Quervains syndrome

FEATURES

- > Plastic structure for excellent immobilization
- > The lightness of the brace and the soft padding grant optimal comfort
- > Adjustable hook and loop closures for easy application

SIZE	WRIST CIRCUMFERENCE
S/M	13 cm - 19 cm
L/XL	19 cm - 27 cm

THUMB WRIST SPLINT**5325-left****THUMB WRIST SPLINT****5326-right**



COOLMAX®

COOLMAX® is a breathable fabric that leaves evaporate moisture from skin making it fresh and dry and ensures maximum comfort to the patient. Coolmax® is a trademark of INVISTA.

INDICATIONS

- > Immobilization of the wrist
- > Sprains
- > Carpal tunnel syndrome
- > Tendinitis
- > Arthritis
- > Contractures

FEATURES

- > In COOLMAX® fabric
- > The neoprene structure allows it to be applied even in the event of a wrist injury
- > Adjustable hook and loop closure to facilitate adjustment

SIZE	WRIST CIRCUMFERENCE
ONE SIZE	up to 28 cm

COOLMAX® WRIST SUPPORT 53300



COOLMAX® DUAL PAD ELBOW STRAP 5345



INDICATIONS

- > Compression of tendon insertions at the elbow
- > Resolution of pain caused by sports injuries and by elbow overexertion (tennis elbow)
- > Acute and chronic epicondylitis

FEATURES

- > COOLMAX® fabric
- > Anatomically shaped neoprene
- > Silicone pads to increase pressure
- > Adjustable hook and loop closure for compression control
- > Elastic finishing edge
- > Ambidextrous

SIZE	ELBOW CIRCUMFERENCE
S	22 cm – 27 cm
M	24 cm – 29 cm
L	26 cm – 31 cm
XL	29 cm – 34 cm

INDICATIONS

- > Finger immobilization and protection
- > Fractures
- > Sprains
- > Dislocations

FEATURES

- > Plastic shell to provide superior protection and immobilization to injured finger
- > Hook and loop closures for an optimal stability
- > Lightweight and of easy application
- > Elastic finishing edge

SIZE	LENGTH OF THE SPLINT
S	5,7 cm
M	6,8 cm
L	8,3 cm
XL	9,3 cm

FINGER SPLINT 5319



SUPER ELASTIC KNEE SUPPORT

5701

INDICATIONS

- > Support and compression for weak or overworked knee
- > Patellar instability
- > Prevention in sports practice

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm



SUPER ELASTIC KNEE SUPPORT with patella opening

5702



INDICATIONS

- > Support and compression for weak or overworked knee
- > Patellar syndrome
- > Patellar tendinopathy
- > Prevention in sports practice

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > Patella opening
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm

INDICATIONS

- > Compression and patellar instability
- > Stabilization of the knee joints
- > Reduction of pain and swelling
- > Consequences of patella dislocation

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > Sides steel splints to provide additional stability
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm

SUPER ELASTIC KNEE BRACE

5703



SUPER ELASTIC KNEE BRACE with patella opening

5704



INDICATIONS

- > Compression and patellar instability
- > Stabilization of the knee joints
- > Reduction of pain and swelling
- > Consequences of patella dislocation

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > With patella opening
- > Sides steel splints to provide additional stability
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm

KNEE IMMOBILIZER 3-panels

20"

57150

INDICATIONS

- > Post-operative treatment of knee ligament reconstructions
- > Post-traumatic phase knee sprains or patella fracture
- > In emergency and first aid situations
- > To replace the casting immobilization

KNEE IMMOBILIZER 3-panels

24"

57160

FEATURES

- > Fabric with soft foam padding
- > Mouldable rigid side and back splints
- > Three adjustable panels
- > Patellar opening
- > Adjustable fastening straps with hook and loop closure
- > Ambidextrous

SIZE

ONE SIZE



ROM KNEE BRACE

19"

57110

ROM KNEE BRACE

22"

57120

INDICATIONS


- > Patella fractures
- > Post traumatic knee sprains
- > Post-operative treatment of knee ligament reconstructions
- > Early removal of casting

FEATURES

- > Aluminium hinges
- > Complete range of motion control for flexion and extension
- > Hook and loop closure straps
- > Soft polyurethane foam padding
- > Ultralight
- > Ambidextrous

SIZE	TIGHT CIRCUMFERENCE	CALF CIRCUMFERENCE
ONE SIZE	up to 70 cm	up to 48 cm





COOLPLUS®

COOLPLUS® is a revolutionary breathable fabric which, thanks to its fiber structure, leaves the skin fresh and dry.
Copyright © CHUNG SHING TEXTILE CO., LTD.



COOLPLUS® MESH HINGED KNEE BRACE

5717

INDICATIONS


- > Instability of the medial and lateral ligaments
- > Prevention and treatment of knee sprain injuries
- > Arthritic instability
- > Senile knee pathology
- > Prevention in sports activity

FEATURES

- > COOLPLUS® breathable fabric
- > Patellar opening with stabilizer
- > Front opening
- > Adjustable hook and loop closures above and below the knee to offer greater stability
- > Aluminum side splints to improve immobilization
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	33 cm – 38 cm
M	38 cm – 43 cm
L	43 cm – 48 cm
XL	48 cm – 53 cm
XXL	53 cm – 58 cm

BREATHABLE



COOLPLUS®

POST-OP ROM KNEE BRACE

13"

5724

COOLPLUS®

POST-OP ROM KNEE BRACE

17"

5725

INDICATIONS




- > Severe knee injuries
- > Post surgery for anterior and posterior cruciate reconstruction
- > Medial and lateral collateral ligament repair
- > Meniscal revisions

FEATURES

- > Made of COOLPLUS® breathable fabric
- > Splints with bilateral ROM mechanism
- > Patented ROM mechanism
- > Flexion-extension stop kit
- > System with extra extension stops
- > Patellar opening with stabilizer
- > Hook and loop adjustable closure straps
- > Ultralight
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	33 cm – 38 cm
M	38 cm – 43 cm
L	43 cm – 48 cm
XL	48 cm – 53 cm
XXL	53 cm – 58 cm

COOLPLUS® is a revolutionary reathable fabric which, thanks to its hair structure, leaves the skin fresh and dry.





COOLMAX®

COOLMAX® is a breathable fabric that leaves evaporate moisture from skin making it fresh and dry and ensures maximum comfort to the patient. Coolmax® is a trademark of INVISTA.

COOLMAX® SPORT KNEE SUPPORT 57200



INDICATIONS

- > Patellar instability
- > Consequences of patella dislocation
- > Prevention in sports practice

FEATURES

- > COOLMAX® breathable fabric
- > Adjustable hook and loop closures to offer greater stability and a comfortable fit
- > Patellar opening with stabilizer
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
ONE SIZE	up to 55 cm

INDICATIONS

- > Patellar syndrome
- > Pre-post accident patella blocking
- > Overloading of the sportsman's patellar tendon
- > Inflammations
- > Tendinopathies

FEATURES

- > COOLMAX® breathable fabric
- > Anatomically shaped orthopedic neoprene
- > With extra pad to provide patella compression
- > Adjustable hook and loop closures to offer greater stability
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
ONE SIZE	up to 55 cm

COOLMAX® KNEE SUPPORT 57300



COOLMAX® PATELLA STRAP with extra pad 57900



INDICATIONS

- > Patellar syndrome
- > Pre-post accident patella blocking
- > Overloading of the sportsman's patellar tendon
- > Inflammations
- > Tendinopathies

FEATURES

- > COOLMAX® breathable fabric
- > Anatomically shaped orthopedic neoprene
- > With extra pad to provide patella compression
- > Adjustable hook and loop closures to offer greater stability
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
ONE SIZE	up to 60 cm

ANKLE SUPPORT

5901

INDICATIONS

- > Instability of the ankle
- > Sprains
- > Edema
- > Prevention of accidental trauma to ankle in sport
- > Useful in re-start of sporting activity

FEATURES

- > Super elastic fabric
- > Heel opening
- > Easy to wear
- > Washable
- > Elastic finishing edge
- > Ambidextrous

SIZE	ANKLE CIRCUMFERENCE
S	18 cm – 23 cm
M	21 cm – 25 cm
L	24 cm – 29 cm
XL	29 cm – 33 cm



COOLMAX® ANKLE STRAP

59200

INDICATIONS

- > From mild to moderate ankle sprains
- > Preventive bandaging

FEATURES

- > COOLMAX® breathable fabric
- > Neoprene bandaging material can be easily applied to swollen or injured ankle
- > Hook and loop closures for a comfortable and easy fit
- > One size fits all

SIZE
ONE SIZE



AIR ANKLE STIRRUP BRACE

10"

59180

AIR ANKLE STIRRUP BRACE

9"

59190

INDICATIONS

- > Prevention and treatment of ankle sprains
- > Tendinitis
- > Prevention of accidental ankle trauma in sporting activity
- > Arthrosis of the ankle joint

FEATURES

- > Structure in composite plastic material, resistant and anatomic
- > Anatomical lateral shells for ankle immobilization
- > Side air cushions to aid circulation and reduce swelling
- > Hook and loop closures for more comfort
- > Ambidextrous

SIZE

ONE SIZE

ANKLE WALKER

17"

5907

INDICATIONS

- > Acute ligament injuries
- > Rupture of the Achilles tendon
- > Malleolar fractures
- > Post-operative treatment of leg and foot fractures
- > Substitute for casting

FEATURES

- > Durable light material construction
- > Wide and soft support base to increase comfort and stability
- > EVA foot bed padding to ensure flexibility
- > The tilting sole helps promote natural stepping and reduce plantar pressure
- > Hook and loop closure straps
- > Ambidextrous

SIZE	MAX FOOT LENGHT	MAX SHOE SIZE
S	24 cm	36 - 37
M	28 cm	38 - 40
L	32 cm	41 - 44
XL	36 cm	45 - 48



ROM ANKLE WALKER

17"

5909



INDICATIONS

- > Acute ligament injuries
- > Rupture of the Achilles tendon
- > Malleolar fractures
- > Post-operative treatment of leg and foot fractures
- > Substitute for casting

FEATURES

- > Durable light material construction
- > ROM mechanism to control range of motion
- > ROM settings between 0° and 45° plantar and dorsal flexion (increments of 15°)
- > Wide and soft support base to increase comfort and stability
- > EVA foot bed padding to ensure flexibility
- > Hook and loop closure straps
- > Ambidextrous

SIZE	MAX FOOT LENGTH	MAX SHOE SIZE
S	24 cm	36 - 37
M	28 cm	38 - 40
L	32 cm	41 - 44
XL	36 cm	45 - 48

INDICATIONS

- > Foot and ankle fractures
- > Joints sprains
- > Acute knee distorsions
- > Post-operative treatment
- > Substitute for casting

FEATURES

- > Rigid construction in durable and lightweight composite material to ensure optimal support and immobilization to the foot and the ankle
- > Wide and soft support base to increase comfort and stability
- > Soft, breathable PU padding for greater comfort
- > The tilting sole helps the natural gait
- > Non-slip and wear-resistant sole for greater safety
- > Ambidextrous

SIZE	BASE LENGTH	BASE WIDTH	HEIGHT	SHOE SIZE
S	16,3 cm	7,8 cm	17,5 cm	22-26
M	20,1 cm	9,1 cm	22,5 cm	27-31
L	22,5 cm	10,4 cm	28,7 cm	32-35

PEDIATRIC WALKER

5927

