



www.medi4care.com

CERTIFICATIONS

The Notified Body selected to certify the Quality Management System of our Company is TÜV Rheinland Italy.

Our certifications represent an important commitment to the Health Care market.



ISO 9001

The certification of the Quality Management System according to ISO 9001-2015 and ISO 13485-2016 standards gives to the Medical Device Manufacturers the warranty that they are producing and marketing their products in accordance to the specific requirements of the market and of the reference laws.



ISO 13485

ISO 13485 is a standard that represents affinity with the concepts defined in the ISO 9001 Norm, increasing its specific value for the medical sector by requiring particular attention to relevant requirements for the safety of Medical Devices placed on the market.



WHO WE ARE

We are a modern and dynamic Company born in 2011 from the entrepreneurial ability and thirty years experience of our Management.

The search of innovative products and new **solutions for the Health Care market** and users is the key element of Medi4care's development and innovation process.

On a daily basis, we are committed to shaping expertise through a wide and accurate range of Medical Devices for Mobility, Home Care, Orthopedics, Rehabilitation, Diagnostics and Wound Care.

OUR MISSION

Improving people's life

Our choices are led by the ideal that is the very essence of the Company "Improving people's life": **to improve the quality of life of people and caregivers** is for us a must.

We are continuously engaged to **the Design, Research and Development** of Medical Devices, offering a choice **of innovative and quality solutions** in order to guarantee a keen service for patients and Health Care professionals.

OUR VALUES

Reliability, Passion, Commitment.

Reliability, professionalism, passion but above all commitment: commitment to corporate social responsibility to always ensure greater attention to the needs of people, community and environment.

Another **value to be shared**, which is an integral part of our way of operating, is hence Sustainability and with this in mind we have adopted Environmentally Friendly models of life: in 2012 a Photovoltaic system was put into operation to cover Medi4Care's energy consumption, to produce energy from renewable sources and to reduce the emission of CO₂ into the atmosphere.

INDEX

MOBILITY

REHALIFE® Wheelchairs	pag. 6
REHALIFE® ALU Wheelchairs	pag. 10
Accessories and spare parts for wheelchairs	pag. 11

HOME CARE

Lift chairs	pag. 14
Antidecubitus	pag. 21
Walking aids	pag. 26
Commode chairs	pag. 33
Bath aids	pag. 36
Daily life aids	pag. 40

BRACES & SUPPORTS - CONWELL

Neck	pag. 45
Shoulder	pag. 47
Trunk	pag. 51
Wrist - Elbow	pag. 55
Knee	pag. 59
Ankle	pag. 65

DIAGNOSTICS

Thermometers	pag. 70
Blood pressure monitors	pag. 71
Pulse oximeters	pag. 73

WOUND CARE

MEDI4SKIN® PU film dressings	pag. 76
MEDI4AKT® HD Hydrocolloids dressings	pag. 82
MEDI4AKT® SF PU Foam and Silicon dressings	pag. 86
MEDI4AKT® F Tracheo - Tracheotomy dressings	pag. 90

FREE TO MOVE



MOBILITY

REHALIFE®

A complete range of standard and light self-propelling **Wheelchairs**,
spare parts and **accessories**, for the disabled person with walking difficulties.

REHALIFE®

Foldable single crossbar steel wheelchair

FEATURES

- > Powder coated steel frame
- > ø 60 cm solid PU rear wheels
- > 20 × 5 cm PU front castors
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Blue color

 SEAT WIDTH	 TOTAL WIDTH	 DIAMETER BACK WHEEL	 DIAMETER FRONT WHEEL	 SEAT HEIGHT	 TOTAL HEIGHT	 WEIGHT	 MAX WEIGHT
41 cm	58 cm	60 cm	20 cm	50 cm	91 cm	18 kg	120 kg
43 cm	60 cm	60 cm	20 cm	50 cm	91 cm	18 kg	120 kg
46 cm	63 cm	60 cm	20 cm	50 cm	91 cm	18,5 kg	120 kg
48 cm	65 cm	60 cm	20 cm	50 cm	91 cm	18,5 kg	120 kg

Ref. 400800
Ref. 400801
Ref. 400802
Ref. 400803



Ref. 400840
Ref. 400841

REHALIFE® FIT

Narrow width foldable single crossbar steel wheelchair

FEATURES

- > Narrow width powder coated steel frame
- > ø 50 cm solid PU rear wheels
- > 20 × 5 cm PU front castors
- > Height adjustable and swing away footrests
- > Removable desk armrests
- > Silver color

 SEAT WIDTH	 TOTAL WIDTH	 DIAMETER BACK WHEEL	 DIAMETER FRONT WHEEL	 SEAT HEIGHT	 TOTAL HEIGHT	 WEIGHT	 MAX WEIGHT
41 cm	52 cm	51 cm	20 cm	50 cm	92 cm	18 kg	120 kg
43 cm	54 cm	51 cm	20 cm	50 cm	92 cm	18 Kg	120 kg

REHALIFE® DOUBLE**Foldable double crossbar steel wheelchair****FEATURES**

- > Powder coated steel frame
- > ø 60 cm solid PU rear wheels
- > 20 × 5 cm PU front castors
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Blue color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
41 cm	58 cm	60 cm	20 cm	50 cm	91 cm	18,5 kg	120 kg
43 cm	60 cm	60 cm	20 cm	50 cm	91 cm	18,5 kg	120 kg
46 cm	63 cm	60 cm	20 cm	50 cm	91 cm	19,0 kg	120 kg
48 cm	65 cm	60 cm	20 cm	50 cm	91 cm	19,0 kg	120 kg

Ref. 400810
Ref. 400811
Ref. 400812
Ref. 400813



Ref. 400815
Ref. 400816

REHALIFE® DOUBLE XL**Bariatric foldable double crossbar steel wheelchair****FEATURES**

- > Powder coated heavy duty steel frame
- > ø 60 cm solid PU rear wheels
- > 20 × 5 cm PU front castors
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Dark silver color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
55 cm	73 cm	60 cm	20 cm	50 cm	91 cm	20,5 Kg	180 kg
60 cm	78 cm	60 cm	20 cm	50 cm	91 cm	20,5 Kg	180 kg



Ref. 400830
Ref. 400831
Ref. 400832
Ref. 400833

REHALIFE® PLUS**Foldable double crossbar steel wheelchair****FEATURES**

- > Powder coated steel frame
- > ø 60 cm solid PU rear wheels with quick release
- > 20 × 5 cm PU front castors
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Safety wheels
- > Blue color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
41 cm	58 cm	60 cm	20 cm	50 cm	91 cm	19 Kg	120 Kg
43 cm	60 cm	60 cm	20 cm	50 cm	91 cm	19 Kg	120 Kg
46 cm	63 cm	60 cm	20 cm	50 cm	91 cm	19,5 kg	120 Kg
48 cm	65 cm	60 cm	20 cm	50 cm	91 cm	19,5 Kg	120 Kg

REHALIFE® COMFORT

Foldable single crossbar steel wheelchair with WC


FEATURES

- > Powder coated steel frame
- > Folding backrest
- > 45 cm width seat with hygiene recess
- > ø 60 cm solid PU rear wheels with quick release
- > 20 x 2,5 cm solid PVC front castors
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Blue color

							
SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
45 cm	63 cm	60 cm	20 cm	54 cm	97 cm	20 kg	125 kg



Ref. 400845



Ref. 400850
Ref. 400851
Ref. 400852
Ref. 400853

REHALIFE® TRANSIT

Transit foldable single crossbar steel wheelchair

FEATURES

- > Powder coated steel frame
- > ø 30 cm solid PU rear wheels
- > 20 x 5 cm PU front castors
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Blue color

							
SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
41 cm	58 cm	30 cm	20 cm	50 cm	91 cm	15 kg	110 kg
43 cm	60 cm	30 cm	20 cm	50 cm	91 cm	15 kg	110 kg
46 cm	63 cm	30 cm	20 cm	50 cm	91 cm	15,5 kg	110 kg
48 cm	65 cm	30 cm	20 cm	50 cm	91 cm	15,5 kg	110 kg

REHALIFE® DOUBLE TRANSIT

Transit foldable double crossbar steel wheelchair

FEATURES

- > Powder coated steel frame
- > ø 30 cm solid PU rear wheels
- > 20 x 5 cm PU front castors
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Blue color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
41 cm	58 cm	30 cm	20 cm	50 cm	91 cm	16,5 kg	120 kg
43 cm	60 cm	30 cm	20 cm	50 cm	91 cm	16,5 kg	120 kg
46 cm	63 cm	30 cm	20 cm	50 cm	91 cm	17 kg	120 kg
48 cm	65 cm	30 cm	20 cm	50 cm	91 cm	17 kg	120 kg

Ref. 400855
Ref. 400856
Ref.400857
Ref.400858



Ref. 400821
Ref. 400822

REHALIFE® EVOLUTION

Convertible foldable single crossbar steel wheelchair

FEATURES

- > Powder coated steel frame
- > 2 positions frame to convert the wheelchair from self-propelling into a transit one
- > ø 60 cm solid PU rear wheels with quick release
- > Equipped with ø 30 cm solid PU rear wheels with quick release for the conversion into a transit wheelchair
- > 20 x 5 cm PU front castors
- > Safety wheels
- > Height adjustable and swing away footrests
- > Flip-up desk armrests
- > Blue color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
43 cm	62 cm	30/60 cm	20 cm	50 cm	88 cm	17 kg	110 kg
46 cm	65 cm	30/60 cm	20 cm	50 cm	88 cm	17 Kg	110 kg

REHALIFE® ALU ONE

Lightweight foldable double crossbar aluminium wheelchair

FEATURES

- > Liquid painting aluminium frame
- > ø 60 cm solid PU rear wheels with quick release
- > 20 x 3,2 cm height adjustable PU front castors
- > Height adjustable and detachable footrests
- > Flip-up PU desk armrests
- > Safety wheels
- > Light blue color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
41 cm	60 cm	60 cm	20 cm	50 cm	92 cm	14 kg	120 kg
43 cm	62 cm	60 cm	20 cm	50 cm	92 cm	14 kg	120 kg
46 cm	65 cm	60 cm	20 cm	50 cm	92 cm	14 kg	120 kg
48 cm	67 cm	60 cm	20 cm	50 cm	92 cm	14 kg	120 kg

Ref. 400900
Ref. 400901
Ref. 400902
Ref. 400903



REHALIFE® ALU SPRINT

Lightweight foldable single crossbar aluminium wheelchair

FEATURES

- > Liquid painting aluminium frame
- > Half foldable backrest
- > Comfort seat and backrest cushion
- > ø 60 cm solid PU rear wheels with quick release
- > 17,5 x 3,2 cm solid PU front castors
- > Height adjustable and detachable footrests with heel straps
- > Height adjustable PU armrests
- > Safety wheels
- > White color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
40 cm	61 cm	60 cm	17,5 cm	54 cm	95 cm	14 kg	125 kg
45 cm	66 cm	60 cm	17,5 cm	54 cm	95 cm	14 kg	125 kg
50 cm	71 cm	60 cm	17,5 cm	54 cm	95 cm	14 kg	125 kg

Ref. 400920
Ref. 400921
Ref. 400922



REHALIFE® ALU PLUS

Lightweight foldable double crossbar aluminium wheelchair

FEATURES

- > Liquid painting aluminium frame
- > ø 60 cm solid PU rear wheels with quick release
- > 20 x 3,2 cm height adjustable PU front castors
- > Height adjustable and detachable footrests
- > Flip-up and height adjustable PU desk armrests
- > Height adjustable safety wheels
- > Light blue color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
41 cm	60 cm	60 cm	20 cm	50 cm	92 cm	15 kg	120 kg
43 cm	62 cm	60 cm	20 cm	50 cm	92 cm	15 kg	120 kg
46 cm	65 cm	60 cm	20 cm	50 cm	92 cm	15 kg	120 kg
48 cm	67 cm	60 cm	20 cm	50 cm	92 cm	15 kg	120 kg

Ref. 400910
Ref. 400911
Ref. 400912
Ref. 400913



REHALIFE® ALU TRANSIT

Lightweight foldable double crossbar aluminium transit wheelchair

FEATURES

- > Liquid painting aluminium frame
- > ø 30 cm solid PU rear wheels
- > 20 x 3,2 cm height adjustable PU front castors
- > Height adjustable and detachable footrests
- > Flip-up PU desk armrests
- > Light blue color

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
41 cm	58 cm	30 cm	20 cm	50 cm	92 cm	13 kg	120 kg
43 cm	60 cm	30 cm	20 cm	50 cm	92 cm	13 kg	120 kg
46 cm	63 cm	30 cm	20 cm	50 cm	92 cm	13 kg	120 kg
48 cm	65 cm	30 cm	20 cm	50 cm	92 cm	13 kg	120 kg

Ref.400940

Ref.400941

Ref.400942

Ref.400943



Ref. 400860

REHALIFE® EASY

Travel foldable aluminium wheelchair

FEATURES

- > Lightweight aluminium frame
- > 5" PVC front castors and 8" PVC rear wheels
- > Flip-up EVA padded armrests
- > Flip-up footrest
- > Lever brake
- > Safety belt with click buckle
- > Transport bag included

SEAT WIDTH	TOTAL WIDTH	DIAMETER BACK WHEEL	DIAMETER FRONT WHEEL	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
30 cm	46 cm	20 cm	12,5 cm	44 cm	95 cm	7 kg	100 kg

ACCESSORIES AND SPARE PARTS

Ref. 409800 60 cm Pneumatic for all the models

Ref. 409801 30 cm Q.R. PU wheels for Rehalife Evolution pair

Ref. 409802 60 cm Pneumatic back wheel for all models

Ref. 409803 60 cm PU back wheel for all models

Ref. 409810 20x3,2 cm PU front wheel for all models of Rehalife Alu

Ref. 409820 20x5 cm PU front wheel for all models of Rehalife

Ref. 409836 Safety wheels for Rehalife Alu One and Rehalife Alu Transit (pair)

Ref. 409837 Safety wheels for Rehalife Alu Plus (pair)

Ref. 409838 IV Stand for wheelchairs

Ref. 409839 Safety wheels for Rehalife (pair)

Ref. 409840 Flip-up desk armrests for all models Rehalife (pair)

Ref. 409845 Flip-up long armrests for all models Rehalife (pair)

Ref. 409850 Flip-up desk armrests for Rehalife Alu One and Alu Transit (pair)

Ref. 409851 Height adjustable armrests for Rehalife Alu Plus (pair)

ACCESSORIES AND SPARE PARTS

Ref. 409860 Shaped desk for wheelchair
Ref. 409861 Wheelchair linear desk
Ref. 409862 Elevating footrests for all models Rehalife (pair)
Ref. 409863 Elevating footrests for all models Rehalife Alu (pair)
Ref. 409870 Swing away trunk supports
Ref. 409874 Abdominal safety belt
Ref. 409875 Abdominal safety belt with braces
Ref. 409877 Abdominal safety belt with pelvic straps and braces

ABDOMINAL SAFETY BELT

FEATURES

- > Breathable TRI-TEX fabric
- > PU foam padding
- > Easy open hook and loop closure
- > Max circumference 130 cm

Ref.409874



ABDOMINAL SAFETY BELT WITH BRACES

FEATURES

- > Breathable TRI-TEX fabric
- > PU foam padding
- > 1 ventral belt and 2 adjustable braces
- > Easy open hook and loop closure
- > Max circumference 130 cm

Ref.409875



ABDOMINAL SAFETY BELT WITH PELVIC STRAPS AND BRACES

FEATURES

- > Breathable TRI-TEX fabric
- > PU foam padding
- > 1 ventral belt, 2 braces and 2 pelvic straps
- > Easy open hook and loop closure
- > Max circumference 130 cm

Ref.409877



IV STAND FOR WHEELCHAIR



Ref. 409838

WHEELCHAIR LINEAR DESK



Ref. 409861

SAFETY WHEELS



Ref. 409836

ELEVATING FOOTRESTS FOR REHALIFE WHEELCHAIR



Ref. 409862

SWING AWAY TRUNK SUPPORTS



Ref. 409870

BY YOUR SIDE **ANYWHERE**



HOME CARE

Home Care aids for disabled and elderly patients with mobility difficulties:

lift chairs, commode chairs, walkers, rollators, crutches, three legs cane, quad cane, canes, supports, IV poles, tables, toilet raisers, shower and bath seats, rehabilitation, devices for the prevention of decubitus, etc.

DUBLINO

Modern design **lift chair with one and two motors**. It is the ideal chair to be adapted to any type of furniture, even to the smallest spaces, to have **maximum relax and stand up without any effort**.



GREY FABRIC
Ref.401000 1 motor
Ref.401010 2 motors



LIGHT BLUE FABRIC
Ref.401001 1 motor
Ref.401011 2 motors



BORDEAUX FABRIC
Ref.401002 1 motor
Ref.401012 2 motors


Class I Medical Device

FEATURES

- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with greater padding for neck and head support
- > Rounded armrests with wide support base
- > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The 1-motor version reclines the backrest and lifts the footrest in a single movement
- > The 2-motors version allows you to separately adjust the inclination of the backrest and the movement of the footrest
- > Maximum backrest inclination 135 °
- > N. 2 rear wheels
- > Handy side pocket for magazines and remote control
- > Gran Foulard in the same fabric and color as the armchair to protect the upholstery

TECHNICAL FEATURES

1-motor version - Electric adjustment in 3 positions



2-motors version - Electric adjustment in 5 positions



Dimensions in cm									
DUBLINO 1 motor	73	106	91	166	49	47	54	35	130
DUBLINO 2 motors	73	106	91	166	49	47	54	40	130

MADRID

Madrid is the **lift chair** with one and two motors which, thanks to its **slim** line, takes up space in width. However, its design does not compromise the comfort guaranteed by the **soft padding** of the backrest and the standard depth of the seat, completely in line with that of other armchairs. Combining **classic and modern** style, Madrid is available in four colours to meet all furnishing needs and tastes.



BORDEAUX MICROFIBER
Ref.401050 1 motor
Ref.401060 2 motors



CORAL MICROFIBER
Ref. 401051 1 motor
Ref. 401061 2 motors



BLUE MICROFIBER
Ref. 401052 1 motor
Ref. 401062 2 motors



OLIVE GREEN MICROFIBER
Ref. 401053 1 motore
Ref. 401063 2 motori

Class I Medical Device

FEATURES

- > Steel and wood frame
 - > Seat padding in non-deformable Polyurethane foam
 - > Ergonomic backrest with extra padding to ensure maximum comfort for people of all stature
 - > Seagull-wing armrests with a wide support base
 - > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The 1-motor version reclines the backrest and lifts the footrest in a single movement
 - > The 2-motor version allows you to separately adjust the inclination of the backrest and the movement of the footrest
 - > Maximum backrest inclination 135 °
 - > N. 2 rear wheels
 - > Handy side pocket for magazines and remote control
 - > Gran Foulard in the same fabric and color as the armchair to protect the upholstery

TECHNICAL FEATURES

1-motor version - Electric adjustment in 3 positions



2-motors version - Electric adjustment in 5 positions



Dimensions in cm									
MADRID 1 motor	72	100	90	165	50	47	50	35	130
MADRID 2 motors	72	100	90	165	50	47	50	40	130

VIENNA

Lift chair with two motors, it has a soft and welcoming high backrest and rest legs joined to the seat, thus ensuring maximum comfort. This electric lift chair is intended both for people with disabilities and for those who simply want to enjoy maximum comfort and relax at home. It has a **Roller System**, wheels system to facilitate the movement of the chair with person on it.



SAND FABRIC
Ref. 401080



ROYAL BLUE FABRIC
Ref. 401081



COFFEE FABRIC
Ref. 401082


Class I Medical Device

FEATURES

- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with extra padding to ensure the maximum comfort for people of all sizes
- > Rounded armrests with large support base
- > "Stand up" function which allows lifting of the person. With this device the person from the sitting position is accompanied up to an almost erect position, facilitating standing
- > The chair is only available in the version with 2 motors: one allows you to adjust the footrest and the other the inclination of the backrest up to 170° reaching the lying position
- > Roller System mechanism for easy movement
- > Handy side pocket for magazine and remote control
- > Gran Foulard in the same fabric and color as the chair to protect the upholstery

TECHNICAL FEATURES

2-motors - Electric adjustment in 6 positions



Dimensions in cm									
VIENNA 2 motors	82	106	92	165	56	49	52	53	170

GINEVRA

Ginevra **lift chair with one and two motors**, with a classic and elegant line that makes it particularly suitable in any type of environment and pleasant for people of all ages. **The extra padding of the backrest** and the **headrest flaps** ensure maximum comfort in all positions. The lift chair has a **Roller System**, a wheel system that allows to move the chair without any effort even with the person seated on it.



COFFEE MICROFIBER
Ref.401020 1 motor
Ref.401030 2 motors



BEIGE MICROFIBER
Ref.401021 1 motor
Ref.401031 2 motors



BLUE MICROFIBER
Ref.401022 1 motor
Ref.401032 2 motors

Class I Medical Device

FEATURES


- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with extra padding to ensure maximum comfort for people of all stature
- > Rounded armrests with a wide support base
- > Headrest on the sides of the backrest for a comfortable support of the head
- > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The 1-motor version lifts the footrest and reclines the backrest up to a maximum of 135 °
- > The 2-motors version allows you to separately adjust the movement of the footrest and the inclination of the backrest up to 170 °, reaching the lying position
- > N. 2 rear wheels (in the version 1 motor only)
- > Handy side pocket for magazines and remote control
- > Gran Foulard in the same fabric and color as the armchair to protect the upholstery
- > Equipped with Roller system mechanism that allows, by pushing with the foot the bar placed behind the chair, to easily move the chair. Once you have entered the mechanism, the lift function is blocked so as to avoid accidental activations

TECHNICAL FEATURES

1-motor version - Electric adjustment in 3 positions



2-motors version - Electric adjustment in 6 positions



Dimensions in cm									
GINEVRA 1 motor	80	107	90	173	54	53	50	46	170
GINEVRA 2 motors	80	107	90	182	54	53	50	51	170

OSLO

Oslo **lift chair** with soft and welcoming outlines is characterized by the presence of many accessories that enrich the design and ease of use. **The movement mechanism** on four independent wheels and the **push handle** allow you to move the chair easily, even with the person seated on it. The **removable headrest flaps** and the **removable armrests** facilitate translation movements from the chair to the bed and vice versa. Equipped with 2 motors; one allows the classic movements "lying down" and "standing up", while the second is intended for the lifting system that allows the user to lift the level of the chair to facilitate the transition from the chair itself to the bed or the wheelchair, eliminating differences in height.



BLACK PU
Ref.401040



IVORY PU
Ref.401041



Class I Medical Device

FEATURES

- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with soft anatomical padding to offer comfort and a correct position
- > Easily removable armrests to facilitate movement and allow the chair to be brought closer to the table
- > Removable headrest flaps for comfortable head support
- > Equipped with the "stand up" function that allows the person to be lifted. With this device, the person is accompanied from a sitting position to an almost erect position, facilitating standing
- > The chair is only available in the 2-motors version: one allows both the backrest to be reclined and the footrest lifted in a single movement, and the "stand up" position; the other allows you to horizontally lift the base of the chair up to a maximum of 18 cm from the floor
- > Maximum backrest inclination 135 °
- > N. 4 swivel wheels of 10 cm
- > Push handle to easily move the person inside the home
- > Handy side pocket for magazines and remote control

TECHNICAL FEATURES

Electric adjustment in 4 positions



Dimensions in cm									
OSLO 2 motors	74	109	87	173	56	53	54	50	170

OSLO VM

The uniqueness of the characteristics of the Oslo lift chair is enriched, in the VM version, by the **vibromassage and heating functions**. The complex integrated vibromassage mechanism performs a gentle massage that **stimulates blood circulation, helps revitalize muscle tone** e gives a pleasant feeling of **well-being and relax**. The quality of the materials, the different programs with which to massage different parts of the body, as well as the possibility of raising and lowering the seat, backrest and footrest, also taking advantage of the area it heats for **soothe tensions or stiffenings**, make the Oslo VM lift chair a choice of election for those who want to combine relax to self-care. A comfortable and practical piece of furniture not only for the elderly or people with mobility difficulties but also for those users who like to relax at the end of the day, combining **comfort** to the **benefits of a massage**.



Class I Medical Device

FEATURES

- > The Oslo VM relaxation massage lift chair, in addition to Oslo features, combines the two aspects that are most popular and that many want to find: a better quality of life and one pleasant sensation of comfort and well-being throughout the whole body
- > Taking advantage of the massaging swing system, through eight vibrating motors, the chair takes care of stimulating several body parts (back, lumbar, buttocks and legs) with a variable power, ideal for all needs
- > At shoulder vibrations facilitate melting of tension in the cervical area allowing relaxation of the neck and back muscles
- > In the lumbar region, thanks to controlled and continuous use of the vibromassage treatment, you can have significant benefits in the treatment of muscle stress
- > In the buttocks the massage promotes blood circulation e decongests the muscle groups; it is particularly suitable for people who lead a sedentary life
- > In the area of the legs, at the height of the calves, the vibromassage promotes blood circulation, and, if applied while the chair is in a horizontal position, facilitates blood flow in the peripheral areas of the body
- > An easy to use multifunction remote control with LCD screen to program the timer, choose from 10 massage programs and 10 different intensity levels or one manual selection that allows the user to customize the massage optimizing the result
- > Thanks to a heating system positioned in the lumbar area, functional for those who, in addition to massage, also appreciate the soothing action of heat, the Oslo VM lift chair relaxes and instills relaxation and pleasure
- > The Polyurethane cover, easily manageable in terms of daily care, gives practicality and elegance to the lift chair

TECHNICAL FEATURES

Electric adjustment in 4 positions									
Dimensions in cm									
OSLO VM 2 motors	74	109	87	173	56	53	54	50	170

ZERO GRAVITY

Lift chair with a classic and elegant taste that make it suitable not only for people with walking difficulties but also for those looking for a relaxing armchair. This chair has **3 motors** that allow independent adjustment of the footrest and inclination of the backrest, as well as to assume the **Trendelenburg position** or **anti-shock position** to get better perfusion of vital organs such as brain, heart and kidneys.



MELANGE BROWN FABRIC
Ref. 401090



MELANGE BLUE FABRIC
Ref. 401091



Class I Medical Device

FEATURES


- > Steel and wood frame
- > Seat padding in non-deformable Polyurethane foam
- > Ergonomic backrest with extra padding to ensure greater support
- > Rounded armrests with large support base
- > "Stand up" function which allows lifting of the person from the sitting position to an almost erect position, facilitating standing
- > Equipped with 3 independent motors
- > N. 4 wheels (the 2 rear ones with brake)
- > Handy side pocket for magazine and remote control
- > Gran Foulard in the same fabric and color as the chair to protect the upholstery

TECHNICAL FEATURES

3-motors - Electric adjustment in 8 positions



Dimensions in cm									
ZERO GRAVITY 3 motors	82	107	92	180	56	49	50	55	170



KIT AIRPUMP WITH PRESSURE ADJUSTEMENT AND BUBBLE MATTRESS MC 2500M

Ref. 400751

FEATURES

- > Simple installation, use and maintenance
- > Pressure regulation depending to the patient weight
- > Alternating cycle operation to prevent bedsores
- > The airpump is fitted with hooks to be positioned at the foot of the bed
- > Mattress consisting of 130 atoxic PVC cells
- > Equipped with end flaps for a better fixing
- > All parts of the mattress are washable
- > Repair kit included

MATTRESS TECHNICAL FEATURES

- > Dimensions (including end flaps) 280 cm (L) x 90 cm (W) x 6.5 cm (H)
- > Beige colour
- > Weight 2 kg
- > Max weight 130 kg

AIRPUMP TECHNICAL FEATURES

- > Power supply 220V / 50 Hz
- > Air output > 4.5 liters per minute
- > Pressure range 45 - 110 mmHg
- > Cycle time 12 minutes
- > Power 7 W
- > 1A Fuse
- > Dimensions 24.5cm (L) x 12cm (W) x 9.3cm (H)
- > Weight 930 g



LIGHT
AND
SIMPLE
TO USE



HANGING
HOOKS



END
FLAPS



130 CELLS



MAX
WEIGHT
130 KG

SPARE PARTS

Bubble mattress 2500M

Ref. 409780



KIT AIRPUMP WITH PRESSURE ADJUSTEMENT AND INTERCHANGEABLE ELEMENTS MATTRESS MC 5000M

Ref. 400752

FEATURES

- > Simple installation, use and maintenance
- > Pressure regulation depending on the patient's weight
- > Alternating cycle operation to prevent bedsores
- > The airpump is equipped with hooks to be positioned at the foot of the bed
- > 17 interchangeable elements in atoxic Nylon / PVC; the tubes can be positioned individually by means of special buttons
- > The first 3 elements are static and they remain inflated to offer support to the cervical area ensuring maximum comfort to the patient
- > Equipped with CPR (Cardio Pulmonary Resuscitation) valve for rapid deflation of the mattress in case of urgent interventions
- > Waterproof and breathable material which does not cause friction and holds all the elements together for better compactness
- > Equipped with end flaps for better fixing
- > All parts of the mattress are washable

MATTRESS TECHNICAL FEATURES

- > Dimensions 190 cm (L) x 88 cm (W) x 11 cm (H)
- > Blue colour
- > Weight 4.7 kg
- > Max weight 140 kg

AIRPUMP TECHNICAL FEATURES

- > Power supply 220V / 50 Hz
- > Air output 5/6 liters per minute
- > Pressure range 40 - 100 mmHg
- > Cycle time 12 minutes
- > Power 7 W
- > 1A Fuse
- > Dimensions 22.4 cm (L) x 10.8 cm (W) x 10.3 cm (H)
- > Weight 1.4 kg



LIGHT
AND
SIMPLE
TO USE



HANGING
HOOKS



END
FLAPS



WATERPROOF
AND
BREATHABLE
MATERIAL




17
INTERCHANGEABLE
ELEMENTS



PILLOW
FUNCTION



MAX
WEIGHT
140 KG



Ref. 400755

3 SECTIONS STATIC PU FOAM MATTRESS

FEATURES

- > The mattress is divided into three sections to fit one or two functions beds giving to the patients the best comfort in all positions
- > The truncated pyramid design ensures air circulation between the mattress and the patient maintaining the ideal environment
- > Density 30 kg/m³

TECHNICAL FEATURES

- > Dimensions: 190 (3 x 63 cm) x 85 x 12,5 cm
- > Max weight 130 kg

ACCESSORIES

Ref. 409790

Ref. 409791

Polyester and Cotton Fabric Cover with zipper 198 x 90 x 15 cm
Waterproof PU Cover with zipper 198 x 90 x 15 cm

MULTI-ZONED STATIC PU FOAM MATTRESS

FEATURES


- > The single section Polyurethane foam mattress is made up of three layers of different densities to ensure differentiated pressure discharge at specific body areas
- > Lateral containment sides density 32 kg/m³
- > Head zone density 29 kg/m³
- > Trunk-pelvis zone density 34 kg/m³
- > Feet perforated zone density 31 kg/m³
- > Equipped with containment cover

TECHNICAL FEATURES

- > Dimensions: 198 x 90 x 15 cm
- > Max weight 130 kg

ACCESSORIES

Polyester and Cotton Fabric Cover with zipper 198 x 90 x 15 cm
Waterproof PU Cover with zipper 198 x 90 x 15 cm



Ref. 400756

Ref. 409790

Ref. 409791


The Air cells cushions Medicare have been designed to protect the patient (wheelchair users, elderly people and people who have to stay in the sitting position for a long period) from the risk of forming pressure sores.




AIR CELLS CUSHION
in PVC and Medical rubber

FEATURES

- > The air cells cushion protects against the risk of pressure sores and tissue injuries
- > The combination of materials, PVC and Medical rubber, allows to maintain the cells soft granting comfort and stability to the patient
- > The micro-interchange cells ensure even weight distribution
- > The air cells cushion is highly conformable to the patient's anatomy whose ischemic points are completely immersed in the cushion favoring regular pressure distribution
- > Cushion cover in breathable elastic fabric that allows air to circulate with every movement of the patient
- > The base of the cushion cover in non-slip material secures the cushion to the chair and helps to increase the stability of the user, guaranteeing thus its therapeutic efficacy
- > Easy and adjustable inflation
- > Lightweight and easy to carry and clean
- > Repair kit and inflation pump included



Pressure level before
the use of the cushion



Pressure level after
the use of the cushion

TECHNICAL FEATURES

REF	HEIGHT	VALVS	DIMENSIONS	COLOUR
400760	6 cm	1	40X40 cm	Grey
400761	6 cm	1	45X40 cm	Grey



LIGHT
AND
SIMPLE
TO USE



MAX
WEIGHT
130 KG

VISCOELASTIC ANATOMIC CUSHION


FEATURES

- > Anatomically shaped cushion in compact viscoelastic polymer which guarantees comfort and stability to the patient and prevents the genesis of pressure sores
- > The elastic PU cover, washable and removable by zipper, is fire retardant and impermeable to liquids
- > Non-slip base cover
- > Equipped with velcro strips for fixing the cushion to the wheelchair and avoiding it slipping

TECHNICAL FEATURES

REF	HEIGHT	DIMENSIONS
400785	6 cm	41X41 cm
400786	6 cm	46X41 cm





Ref. 400775

Ref. 400776

CONTOUR CUSHION WITH SELF-MODELING GEL

FEATURES

- > Composite anti-decubitus cushion with preformed anatomical base integrated with self-modeling gel, indicated for patients who use wheelchairs and long-term patients
- > Made up of viscoelastic "slow memory" foam and fluid gel to reduce the pressure on contact areas most at risk of decubitus
- > Removable PU cover, stretch, washable, with base in non-slip material
- > Equipped with handle to facilitate transport

TECHNICAL FEATURES

- > Dimensions 41x41x7 cm
- > Dimensions 44x44x7 cm

CONTOUR CUSHION WITH SELF-MODELING FLUIDS


FEATURES

- > Composite anti-decubitus cushion with preformed anatomical base integrated with self-modeling fluids, suitable to prevent the formation of pressure sores in patients who use wheelchairs and long-term patients who need a comfortable and effective support
- > Made up of a pre-formed base in compact, high density and high resilience foam, and a layer of viscoelastic fluid gel with a cooling effect
- > Removable PU cover, stretch, washable, with base in non-slip material
- > Equipped with handle to facilitate transport

TECHNICAL FEATURES

Ref. 400780 > Dimensions 41x41x7 cm

Ref. 400781 > Dimensions 44x44x7 cm



Ref. 400770

Ref. 400771

FLUID GEL CUSHION

FEATURES

- > The fluid gel cushion guarantees the optimization of pressures to limit the risks of prolonged capillary occlusion in people with reduced mobility
- > The interior is made up of fluid gel bag enclosed between 2 layers of expanded Polyurethane
- > Equipped with stretch, washable, non-slip removable PU cover

TECHNICAL FEATURES

- > Dimensions 40x40x5 cm
- > Dimensions 45x40x5 cm



Ref. 400765

Ref. 400766

SILICONE COTTON WOOL CUSHION

FEATURES

- > Three-section cushion for the prevention of pressure ulcers
- > The silicone cotton wool padding guarantees softness and breathability to the cushion
- > Covered on both sides with 100% sanforized cotton fabric
- > Equipped with double velcro to secure the cushion to the chair and/or to the wheelchair

TECHNICAL FEATURES

- > Dimensions 40x40x10 cm
- > Dimensions 45x40x10 cm

SILICONE COTTON WOOL CUSHION WITH HOLE

FEATURES

- > Cushion with central hole for the prevention of pressure ulcers in particularly sensitive areas and for prostatic problems
- > The silicone cotton wool padding guarantees softness and breathability to the cushion
- > The 100% sanforized cotton fabric cover allows the skin to stay dry and cool

TECHNICAL FEATURES

- Ref. 400769 > Dimensions 42x42x10 cm



Ref. 400791

SILICONE COTTON WOOL HEEL PROTECTOR

FEATURES

- > The heel protection guarantees to the patient the optimization of the pressures on riskful areas of the body
- > The silicone cotton wool padding reduces the heel pressure and avoids friction and consequent skin irritation
- > The 100% sanforized cotton fabric cover allows the skin to be dry and fresh
- > Equipped with a practical adjustable velcro closure

TECHNICAL FEATURES

- > Dimensions 14x10x3 cm

SILICONE COTTON WOOL ELBOW PROTECTOR

FEATURES

- > The elbow protection guarantees to the patient the optimization of pressures on riskful areas of the body
- > The silicone cotton wool padding reduces the elbow pressure and avoids friction and consequent skin irritation
- > The 100% sanforized cotton fabric cover allows the skin to be dry and fresh
- > Equipped with a practical adjustable velcro closure

TECHNICAL FEATURES

- Ref. 400790 > Dimensions 30x15x11 cm



WOODEN CANE

FEATURES

- > Height 93 cm
- > Non-slip rubber tip
- > Available in brown and black



Ref. 400251	Curved handle wooden cane	ø 22mm	brown
Ref. 400252	Curved handle wooden cane	ø 22mm	black
Ref. 400271	Hook handle wooden cane	ø 25mm	brown
Ref. 400272	Hook handle wooden cane	ø 25mm	black
Ref. 400261	T-handle wooden cane	ø 22mm	brown
Ref. 400262	T-handle wooden cane	ø 22mm	black



FOLDABLE ALUMINIUM CANE

FEATURES

- > Anodized aluminium structure
- > 10 levels of height adjustment
- > Ergonomic wooden handle
- > Non-slip tip

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
70-93 cm	0,30 kg	100 kg

Ref. 400235



ADJUSTABLE ALUMINIUM CANE

FEATURES

- > Anodized aluminium structure
- > 10 levels of height adjustment
- > Ergonomic wooden handle
- > Non-slip tip
- > Available in 4 patterns

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
70-93 cm	0,30 kg	100 kg



Ref. 400240	Bon ton
Ref. 400241	Optical
Ref. 400242	Flowers
Ref. 400243	Fantasy



CRUTCHES

FEATURES

- > Structure in light anodized aluminium
- > 10 levels of height adjustment
- > Forearm support
- > Ergonomic handgrips with reflectors on the front
- > Non-slip rubber tip
- > Available in different colours

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
90-112 cm	0,40 kg	120 kg

Ref. 400201	grey
Ref. 400202	blue
Ref. 400203	red
Ref. 400204	green



CRUTCHES with soft handle

FEATURES

- > Structure in light anodized aluminium
- > 10 levels of height adjustment every 2,5 cm
- > Forearm support
- > Ergonomic, non-slip and soft handgrips with reflectors on the front
- > Anti-noise ring on the sliding tube
- > Non-slip rubber tip
- > Available in grey and black colours

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
92-116 cm	0,45 kg	120 kg

Ref. 400211	gray
Ref. 400212	black





Ref. 400230 **boy**
Ref. 400231 **adult**

UNDERARM CRUTCHES

FEATURES

- > Structure in light anodized aluminium
- > 9 levels of height adjustment
- > Ergonomic forearm support and handle
- > Non-slip rubber tip

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
97-114 cm	0,80kg	110 kg
113-134 cm	0,90kg	110 kg



THREE LEGS CANE

FEATURES

- > Structure in chromed aluminium
- > Height adjustable
- > Ergonomic handle
- > Non-slip rubber tip
- > Reinforced base

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
74-97 cm	0,76 kg	113 kg

Ref. 400301

QUAD CANE

FEATURES

- > Structure in chromed aluminium
- > Height adjustable
- > Ergonomic handle
- > Non-slip rubber tips
- > Reinforced base

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
74-97 cm	0,78 kg	113 kg

Ref. 400351



HD QUAD CANE

FEATURES

- > Wide base structure in anodized aluminium
- > Height adjustable
- > Heavy duty
- > Ergonomic handgrips
- > Non-slip rubber tips

HEIGHT ADJUSTABLE	WEIGHT	MAX WEIGHT
68-91 cm	0,8 kg	136 kg

Ref. 400355 **green**
Ref. 400356 **bordeaux**



Ref. 400601

WALKER reciprocating

FEATURES

- > Aluminium reciprocating walker
- > Foldable by one button
- > Ergonomic handgrips
- > 8 levels height adjustment
- > 4 non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
78-96 cm	47 X 46 cm	2,4 kg	136 kg



WALKER foldable fix

FEATURES

- > Aluminium fix walker
- > Foldable by two buttons
- > Ergonomic handgrips
- > 8 levels height adjustment
- > 4 non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
79-96 cm	48 X 44 cm	2,6 kg	136 kg

Ref. 400501

WALKER double handgrips

FEATURES

- > Aluminium fix walker
- > Foldable by one button
- > 2 height handgrips to ease the patient to stand up
- > Ergonomic handgrips
- > Height adjustable
- > 4 non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
74-87 cm	56 X 52 cm	2,7 kg	136 kg



Ref. 400510



Ref. 400551

WALKER 2 wheels - 2 tips

FEATURES

- > 2 wheels and 2 tips aluminium walker
- > Foldable by one button
- > Ergonomic handgrips
- > 8 levels height adjustment
- > 2 non-slip rubber tips and 2 x 5" wheels

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
79-95 cm	61 X 52 cm	3 kg	136 kg

WALKER
2 castors - 2 wheels with brake

FEATURES

- > 2 x 3" castors and 2 x 3" wheels with brake aluminum walker
- > Foldable by one button
- > Ergonomic handgrips
- > Height adjustable

HEIGHT ADJUSTABLE	BASE	WEIGHT	MAX WEIGHT
83-100 cm	61 X 55 cm	3 kg	110 kg

Ref. 400552



HIGH WALKER
with underarm supports

FEATURES

- > Foldable high walker with height adjustable soft PU underarm supports
- > Steel frame with 4 castors (2 of them with brake)
- > Equipped with padded seat

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
100-140 cm	78 X 55 cm	60 cm	15 kg	120 kg

Ref. 400402



Ref. 400401

HIGH WALKER
with table

FEATURES

- > Foldable high walker with padded table and handles
- > Steel frame with 4 x 100mm castors (2 of them with brake)
- > Height adjustable table and handles
- > Equipped with padded seat, backrest and object holder hook

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
110-135 cm	73 X 64 cm	55 cm	11 kg	136 kg

WALKER
ACCESSORIES

5" CASTOR EXTENSION LEGS

Ref. 409831



3" CASTOR EXTENSION LEGS

Ref. 409830



5" WHEEL EXTENSION LEGS

Ref. 409833



3" WHEEL EXTENSION LEGS

Ref. 409832



PRESSURE BRAKE
EXTENSION LEGS

Ref. 409835



PRESSURE BRAKE AND
3" WHEELS EXTENSION LEGS

Ref. 409834



ROLLATOR foldable 2 wheels

FEATURES

- > Foldable steel frame rollator with 2 x 5" wheels and 2 tips
- > Ergonomic handgrips
- > Height adjustable
- > Padded seat

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
84-94 cm	61 X 55 cm	55-65 cm	6 kg	120 kg

Ref. 400605



Ref. 400606

ROLLATOR SKINNY foldable 2 wheels

FEATURES

- > Foldable aluminium frame rollator
- > 2 x 6" front wheels and 2 back tips
- > Easily foldable for storage and transport
- > Ergonomic handgrips
- > Padded backrest band
- > Height adjustable in 5 different positions
- > Removable shopping bag

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
78-97 cm	60 X 56cm	47-54 cm	5.5 Kg	136 Kg

ROLLATOR foldable 4 wheels with seat

FEATURES

- > Foldable aluminium frame rollator with 8" front castors and back wheels
- > Handle brake on the back wheels
- > Ergonomic handgrips
- > Height adjustable handles
- > Padded seat and backrest band
- > Equipped with basket and cane holder

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
81-92 cm	61 X 61 cm	54 cm	7.5 kg	120 kg

Ref. 400610



ROLLATOR CUBE
foldable and compact

FEATURES

- > Elegant, light, foldable aluminium frame rollator
- > Oval frame completely foldable to grant minimum storage dimensions (63x30x35 cm) and to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Handle brake
- > Height adjustable ergonomic handgrips for extra comfort
- > Equipped with seat, backrest, shopping bag with reflector and cane holder

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
80-93 cm	63X68 cm	54 cm	7,5 kg	136 kg



Ref. 400620



Ref. 400615

ROLLATOR ECO-CUBE
foldable and compact

FEATURES

- > Elegant foldable rollator which combines the lightness of aluminium frame to the refinement of the wood-effect coloring
- > Oval frame completely foldable to grant minimum storage dimensions (90x26,5x33 cm) and to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Handle brake
- > Height adjustable ergonomic handgrips for extra comfort
- > Equipped with seat, backrest, and cane holder
- > Shopping bag with reflector

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
78-93 cm	69X61,5 cm	55,5 cm	5,5 Kg	136 kg





Ref. 400625

ROLLATOR URBAN 2 functions 4 wheels

FEATURES

- > Practical and comfortable aluminium frame rollator convertible into transit wheelchair
- > Double crossbar foldable frame to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Handle brake
- > Height adjustable ergonomic handles
- > Equipped with seat, backrest, shopping bag with reflector and cane holder
- > Flip-up footrests with heel strip

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
84-97 cm	63X75 cm	55.5 cm	9.2 Kg	136 kg

ROLLATOR URBAN AB with PU forearm pads

FEATURES

- > Height adjustable (5 positions) aluminium frame rollator
- > Side-to-side foldable frame to grant minimum storage dimensions (23x61x100 cm) and to ease transport
- > 8" front castors and back wheels with reflectors which make it particularly suitable for outdoor use
- > Height adjustable soft ergonomic forearm pads
- > Handle steel brakes
- > Equipped with fabric seat and cane holder
- > Detachable shopping bag with reflector

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
100-115 cm	70X61 cm	60.5 cm	7.6 Kg	136 Kg

Ref. 400626



SPARE PARTS

Ref. 409825	Wheels for rollator Ref. 400551	Ref. 409915	Rubber tips for Ref. 400301
Ref. 409827	100 mm back wheels with brake for rollator Ref. 400402 (pair)	Ref. 409920	Rubber tips for Ref. 400351
Ref. 409828	100 mm front wheels for rollator Ref. 400402 (pair)	Ref. 409925	Rubber tips for Ref. 400201-202-203-204
Ref. 409900	Rubber tips for Ref. 400551	Ref. 409926	Rubber tips for ø 22mm wooden canes
Ref. 409905	Rubber tips for Ref. 400601	Ref. 409927	Rubber tips for ø 25mm wooden canes
Ref. 409910	Rubber tips for Ref. 400501	Ref. 409930	PU handgrip for Ref. 400501-400551-400601

COMMODE CHAIR WITH
RECLINING BACKREST
self-propelling

FEATURES

- > Steel frame
- > Reclining and removable padded backrest, 56 cm high
- > Padded seat with anatomic opening for toilet
- > Made up of washable and fire retardant material
- > 4 PU wheels, 2 x 24" (60 cm) and 2 x 8" (20 cm)
- > Brake on the rear wheels
- > Padded, flip-up and removable armrests
- > Height adjustable elevating footrests with padded legrests
- > Flip-up and lenght adjustable footrests
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
66 cm	43 cm	50 cm	106 cm	25 kg	130 kg

HEADREST
for commode chairs



Ref. 409730

Ref. 400050



Ref. 400053

COMMODE CHAIR WITH
RECLINING BACKREST
4" wheels

FEATURES

- > Steel frame
- > Reclining and removable padded backrest, 56 cm high
- > Padded seat with anatomic opening for toilet
- > Made up of washable and fire retardant material
- > 4 x 4" swivel wheels, the 2 rear ones with brake
- > Padded, flip-up and removable armrests
- > Height adjustable elevating footrests with padded legrests
- > Flip-up and lenght adjustable footrests
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
56 cm	43 cm	50 cm	106 cm	21 Kg	130 Kg

COMMODE CHAIR WITH RECLINING BACKREST 8" wheels

FEATURES

- > Steel frame
- > Reclining and removable padded backrest, 56 cm high
- > Padded seat with anatomic opening for toilet
- > Made up of washable and fire retardant material
- > 4 x 8" swivel wheels, the 2 front ones with brake
- > Padded, flip-up and removable armrests
- > Height adjustable elevating footrests with padded legrests
- > Flip-up and length adjustable footrests
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
62 cm	43 cm	50 cm	106 cm	22 kg	130 Kg



Ref. 400055



COMMODE CHAIR

FEATURES

- > Height adjustable steel frame
- > Flip-up armrests
- > Padded seat
- > Water-repellent leather seat
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
55 cm	43.5 cm	50-57 cm	95-103 cm	10 kg	136 kg

Ref. 400034

COMMODE CHAIR with wheels

FEATURES

- > Steel frame with 4 x 5" swivel wheels, 2 of which with brake
- > Flip-up armrests
- > Adjustable and removable footrests
- > Padded seat
- > Water-repellent leather seat
- > WC made of easily washable plastic material

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
60 cm	43 cm	50 cm	88 cm	12 kg	120 kg



Ref. 400035



3 IN 1 COMMODE CHAIR

FEATURES

- > 3 uses in 1: toilet raiser, commode chair, toilet seat
- > Steel frame
- > Height adjustable
- > Plastic armrest and round bucket with lid
- > Non-slip rubber tips

HEIGHT ADJUSTABLE	BASE	SEAT HEIGHT	WEIGHT	MAX WEIGHT
57-72 cm	47 X 58 cm	37-52 cm	4,8 Kg	136 kg

Ref. 400030

SHOWER AND COMMODE CHAIR self-propelling

FEATURES

- > Reinforced anti-rust and anti-corrosion aluminium frame
- > Lightweight and versatile: it can be used as a toilet seat
- > PU seat with anatomic opening for toilet to allow the use of the bucket
- > 2 x 24" PU wheels and 2 x 5" wheels with brake
- > Removable PU backrest
- > Flip-up armrests
- > Swiveling, removable footrests with height-adjustable footplates

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
64 cm	46 cm	50 cm	99 cm	18 kg	120 kg

Ref. 400040



SHOWER AND COMMODE CHAIR 5" wheels

FEATURES

- > Reinforced anti-rust and anti-corrosion aluminum frame
- > Lightweight and versatile: it can be used as a toilet seat
- > PU seat with anatomic opening for toilet to allow the use of the bucket
- > 4 x 5" PU wheels with brake
- > Removable PU backrest
- > Flip-up armrests
- > Swiveling, removable footrests with height-adjustable footplates

TOTAL WIDTH	SEAT WIDTH	SEAT HEIGHT	TOTAL HEIGHT	WEIGHT	MAX WEIGHT
55 cm	46 cm	50 cm	99 cm	15 kg	120 kg

Ref. 400045





Ref. 400101

SHOWER CHAIR

FEATURES

- > Aluminium frame
- > Ergonomic seat with holes for water drainage
- > Side handles
- > Height adjustment at 8 different levels
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	SEAT HEIGHT	WEIGHT	MAX WEIGHT
30X50 cm	37-50 cm	1.8 kg	113 Kg

WALL SHOWER SEAT

FEATURES

- > Flip-up aluminium frame to be fixed to the wall
- > Floor legs
- > Ergonomic seat with holes for water drainage
- > Side handles
- > Height adjustable
- > Materials resistant to corrosion

SEAT	SEAT HEIGHT	WEIGHT	MAX WEIGHT
30X50 cm	40-55 cm	4.8 kg	113 Kg



Ref. 400102



Ref. 400110

SHOWER CHAIR WITH BACKREST

FEATURES

- > Aluminium frame
- > Ergonomic seat and backrest with holes for water drainage
- > Side handles
- > Height adjustment at 8 different levels
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	BACK	SEAT HEIGHT	BACK HEIGHT	WEIGHT	MAX WEIGHT
30X50 cm	35X50 cm	37-50 cm	70-87 cm	2.7 kg	120 kg

FOLDABLE SHOWER CHAIR

FEATURES

- > Aluminum foldable frame
- > Ergonomic backrest and seat with holes for water drainage
- > Height adjustable
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	BACK	SEAT HEIGHT	BACK HEIGHT	WEIGHT	MAX WEIGHT
40x41 cm	20x48 cm	46-54 cm	88-96 cm	8,2 Kg	113 Kg

Ref. 400109



Ref. 400115

SWIVEL SHOWER CHAIR

FEATURES

- > Steel frame
- > Ergonomic swivel PVC seat
- > Armrests for more comfort

SEAT	BACK	BASE	WEIGHT	MAX WEIGHT
40X37 cm	40X37 cm	73X52 cm	4,3 kg	120 kg

ROUND SHOWER SEAT

FEATURES

- > Aluminium frame
- > Round seat
- > Height adjustment at 8 different levels
- > Non-slip rubber tips
- > Materials resistant to corrosion

SEAT	SEAT HEIGHT	WEIGHT	MAX WEIGHT
Ø32 cm	38-55 cm	1,6 kg	100 kg

Ref. 400100



BATH BOARD WITH HANDLE AND ROTATING SEAT

FEATURES

- > High resistance PVC seat
- > Non-slip rotating seat
- > Safety handle and adjustable brackets for a better fitting to the bathtub's dimensions

BASE	WEIGHT	MAX WEIGHT
75X36 cm	1,5 kg	120 kg



Ref. 400125



Ref. 400130

SAFETY HANDLE FOR BATHTUB

FEATURES

- > Steel safety handle to facilitate entry and exit from the bathtub
- > Fixing clamp to the bathtub (maximum width 17,8 cm)

WEIGHT	MAX WEIGHT
2,8 kg	120 kg

TRANSFER SEAT FOR BATHTUB

FEATURES

- > Aluminium frame
- > Non-slip PVC seat
- > Backrest with handle
- > Height adjustable
- > Non-slip suction cup tips
- > Materials resistant to corrosion

SEAT	BACK	SEAT HEIGHT	BACK HEIGHT	WEIGHT	MAX WEIGHT
70X41 cm	35X45 cm	46-56 cm	81-91 cm	4,4 kg	120 kg



Ref. 400120

TOILET RAISER

FEATURES

- > Ergonomic seat for more comfort
- > Washable PVC material with high resistance
- > Front and rear opening to ease accessibility for personal hygiene
- > Fixing system with adjustable side brackets



Ref. 400020 - H 10 cm
Ref. 400025 - H 15 cm



Ref. 400010 - H 10 cm
Ref. 400015 - H 15 cm

TOILET RAISER with lid

FEATURES

- > Ergonomic seat for more comfort
- > Washable PVC material with high resistance, with lid
- > Front and rear opening to ease accessibility for personal hygiene
- > Fixing system with adjustable side brackets



Ref. 400005

TOILET RAISER with armrests

FEATURES

- > Ergonomic seat
- > Simple and safe front fixing system
- > It raises the toilet level of 5"
- > Washable PVC material
- > Front opening to ease accessibility for personal hygiene



Ref. 400028

SOFT TOILET RAISER H 10 cm

FEATURES

- > Ergonomic seat for more comfort
- > Washable soft PU material with high resistance
- > Front opening to ease accessibility for personal hygiene
- > Fixing system with adjustable side brackets



IV STAND

FEATURES

- > Steel frame
- > Height adjustable
- > Equipped with 4 hooks
- > Base with 5 castors

HEIGHT ADJUSTABLE	MAX WEIGHT
133-223 cm	2kgxhook

Ref. 400651



Ref. 400701

BACKREST

FEATURES

- > Bed backrest
- > Steel frame
- > Adjustable in 5 different positions
- > Made of breathable nylon with padded headrest
- > Blue color

BASE	HEIGHT ADJUSTABLE	MAX WEIGHT
54x58 cm	30-54 cm	120 kg

EXAMINATION COUCH

FEATURES

- > Steel frame 190(L) x 65(W) x 75(H) cm
- > PU foam padding and washable leather
- > Adjustable in 5 different positions
- > Roll holder
- > Adjustable backrest

MAX WEIGHT
180 kg

Ref. 400660





Ref. 400951

OVERBED TABLE double tilting desk top

FEATURES

- > Steel frame
- > Double tilting desk top surrounded by edges
- > Height adjustable
- > With 4 swivel wheels, 2 of which with brake
- > Dark brown color

SMALL DESK TOP	BIG DESK TOP	HEIGHT ADJUSTABLE
40 X 20 cm	40 X 68 cm	73-125 cm

OVERBED TABLE tilting desk top

FEATURES

- > Steel frame
- > Tilting desk top
- > Height adjustable
- > With 4 swivel wheels, 2 of which with brake
- > Brown color

DESK TOP HEIGHT ADJUSTABLE	HEIGHT ADJUSTABLE
41 X 65 cm	71-110 cm

Ref. 400952



Ref. 400955

TILTING OVERBED TABLE pneumatic adjustment

FEATURES

- > Steel frame
- > Tilting (0° to 330°) desk top, washable, fire retardant, scratch resistant and highly durable
- > Pneumatic height adjustment
- > With 4 swivel wheels, 2 of which with brake
- > Gray color

TABLE TOP HEIGHT ADJUSTABLE	HEIGHT ADJUSTABLE
48 X 78 cm	76-117 cm

PEDAL EXERCISER

FEATURES

- > Steel frame
- > Strength adjustable pedal
- > Triangle anti-slips tips

Ref.400630



Ref.400631

PEDAL EXERCISER digital display

FEATURES

- > Steel frame
- > Digital display for timer, counting number of repetitions, speed and consumed calories
- > Triangle non-slip tips for more stability and safety

SWIVEL STOOL

FEATURES

- > Plastic base
- > Padded PU leather seat
- > Height adjustable

Ref.400668



Ref.400665

STEPPER

FEATURES

- > Steel frame
- > 2 steps
- > Non-slip footrests
- > Non-slip tips

WHEN A **VALID SUPPORT** IS NEEDED



SUPPORTS AND BRACES

The partnership with Conwell Medical Co. Ltd., a company specialized in design and manufacturing of a wide range of support and braces, is the result of a keen research. The synergy between the two companies results into a **complete range of patented and professional solutions** for the body parts at risk of injury.

CONWELL Medical is a fast growing company specialized in the field of orthopedic braces and supports. Founded in Taiwan in 2005, Conwell Medical Co., Ltd. is a manufacturer, distributor and exporter of its patented brand Conwell. Conwell offers a complete range of orthopedic products for the major injury-prone areas of the body.


"Be Well with Conwell" is the mission of the Company. To accomplish that, Conwell rigorously selects the materials and constantly innovate new designs to ensure that the products provide effectiveness and comfort to customers.

Product quality control is a priority.

In 2009, Conwell received the quality standard ISO 9001 and ISO 13485.

**To keep its products to a level
of excellence CONWELL has chosen
COOLMAX® and COOLPLUS®**

MOISTURE MANAGEMENT + BREATHABILITY = TOTAL COMFORT



COOLMAX® fabric


Why choosing COOLMAX® fabric?

For the simple reason of staying dry and comfortable when you start to perspire. Sweating is one of the body's natural mechanisms to stay cool. Engineered with special technology to wick moisture and enhance drying, COOLMAX® fabric helps you stay dry and energized - exactly what you need to feel comfortable not only during the day, but also through your toughest workouts.

Today COOLMAX® fabric appeals to the fashion conscious who prefer a cool, soft and high performance fabric that fits every lifestyle. COOLMAX® fabric is the confident, sophisticated choice for active consumers.

- COOLMAX® fabric pulls moisture away from your skin.
- COOLMAX® fabric absorbs and spreads moisture out across fabric to enhance Moisture Vapor Transmission rate.
- Air moves in to keep your body cool and dry.

COOLMAX® is a trademark of INVISTA.



COOLPLUS® fabric

Why choosing COOLPLUS® fabric?

It keeps you cool and dry.

COOLPLUS® is a new fiber that features water diffusion and wicking properties. These fibers have excellent texture yielded several of micro-slits within single fiber. As a result the fabric does not stick to the skin and shows cool crispy wearing comfortable due to its capillary action.

In addition, the dry touch can be maintained even when the part is still wet.

INDICATIONS

- > Mild sprains of the cervical spine
- > Stiff neck
- > Cervical arthrosis

FEATURES

- > Made of polyurethane foam
- > Covered with stockinette for maximum patient comfort
- > Hook and loop closure for ease of adjustment
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
S	36 cm – 42 cm	7 cm
M	39 cm – 45 cm	7 cm
L	41 cm – 47 cm	9 cm
XL	44 cm – 50 cm	9 cm

**SOFT
CERVICAL
COLLAR****5101****PLASTIC
ADJUSTABLE
COLLAR****5103****INDICATIONS**

- > Support for minor sprains of the cervical spine
- > Stiff neck
- > Discopathies
- > Cervical arthrosis
- > Post-operative
- > Vertiginous crises of cervical origin

FEATURES

- > Semi-rigid polypropylene structure
- > Foam padding with ventilation holes for comfort of the patient
- > Adjustable in height and circumference
- > Hook and loops closure for easy adjustment
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
S	29 cm – 34 cm	6 cm
M	34 cm – 39 cm	7 cm
L	39 cm – 44 cm	9 cm
XL	44 cm – 49 cm	10 cm

INDICATIONS

- > Support for minor sprains of the cervical spine
- > Stiff neck
- > Discopathies
- > Cervical arthrosis
- > Post-operative
- > Vertiginous crises of cervical origin

FEATURES

- > Semi-rigid polypropylene structure
- > Foam padding with ventilation holes for comfort of the patient
- > Adjustable in height and circumference
- > Hook and loops closure for easy adjustment
- > Special design for excellent support
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
S	29 cm – 34 cm	6 cm
M	34 cm – 39 cm	7 cm
L	39 cm – 44 cm	9 cm

**PLASTIC
ADJUSTABLE
COLLAR
with chin
support****5104**

CERVICAL COLLAR with trachea opening h. 4 1/4"

5108

CERVICAL COLLAR with trachea opening h. 5 1/4"

5109



INDICATIONS

- > Support for minor cervical sprains
- > Stiff neck
- > Discopathies
- > Cervical arthrosis
- > Post-operative
- > Vertiginous crises of cervical origin
- > Useful in case of tracheotomy

FEATURES

- > Lightweight and compact bivalve structure
- > Large trachea opening for easy emergency tracheotomies
- > Hook and loop side closures for an easy fit
- > X-ray transparent
- > Available in different heights

SIZE	NECK CIRCUMFERENCE
S	25 cm – 33 cm
M	33 cm – 40 cm
L	40 cm – 48 cm
XL	48+ cm

DELUXE ADJUSTABLE BIVALVE CERVICAL COLLAR

51100

INDICATIONS

- > Cervical trauma
- > Cervical hernias
- > Instability of the cervical tract due to neurological pathologies
- > Cervicoarthrosis
- > Discoarthrosis
- > Myogenic stiff necks
- > Post-operative immobilization

FEATURES

- > Bivalve design made with light and quality materials to fit any morphological profile and improve patient compliance especially in cases of prolonged use
- > Simple height adjustment: this allows the collar to vary the traction on the patient's cervical area
- > The sternal support of the collar increases the stability of the device
- > Wide trachea opening for easy emergency tracheotomies
- > Soft foam padding prevents pressure sores and grants comfort
- > X-ray transparent

SIZE	NECK CIRCUMFERENCE	HEIGHT
ONE SIZE	21 cm – 49 cm	Adjustable



CLAVICLE BRACE
5201

INDICATIONS

- > Immobilization of the clavicle
- > Subluxation of the joint acromio-clavicular

FEATURES

- > Soft foam padding for patient comfort
- > Padded closed cell foam axillary straps covered with breathable stockinette
- > Hook and loop closure for ease of adjustment

SIZE	CHEST CIRCUMFERENCE
XS	44 cm - 54 cm
S	54 cm - 64 cm
M	64 cm - 76 cm
L	76 cm - 89 cm
XL	89 cm - 104 cm



SHOULDER SUPPORT

5202



INDICATIONS

- > Immobilization of the shoulder
- > Sprains
- > Tendinitis
- > Bursitis

FEATURES

- > Soft foam padded shoulder straps for patient comfort
- > Length-adjustable
- > Ambidextrous

SIZE	ARM CIRCUMFERENCE
S/M	23 cm – 32 cm
L/XL	31 cm – 40 cm

DELUXE SHOULDER IMMOBILIZER

52070

INDICATIONS

- > Immobilization of the shoulder
- > Dislocations
- > Subluxations
- > Post-surgery immobilization

FEATURES

- > Soft foam padded bands to ensure comfort to the patient
- > Two over-the-shoulder criss-cross straps to increase joints stability
- > Length-adjustable straps
- > Final stabilization band with hook and loop closure
- > Hook and loop closure closures for an easy fit
- > Ambidextrous

SIZE

ONE SIZE



SHOULDER IMMOBILIZER 5206

INDICATIONS

- > Treatment of dislocations
- > Post-surgery immobilization

FEATURES

- > Elastic waste band construction for compressive support
- > Durable wrist and humeral cuff with contact closure for an easy fit
- > Ambidextrous

SIZE	WAIST CIRCUMFERENCE
S/M	66 cm – 89 cm
L/XL	89 cm – 112 cm



MESH ARM SLING 5204

INDICATIONS

- > Immobilization and support for mild spains and strains
- > Dislocations of the shoulder complex
- > Post fractures

FEATURES

- > Envelope design with thumb support
- > Hook and loop closure for ease of adjustment
- > Ambidextrous

SIZE	ARM LENGHT (from elbow to wrist)
S	28 cm – 33 cm
M	33 cm – 38 cm
L	38 cm – 43 cm
XL	43 cm – 48 cm



DELUXE ARM SLING 5205

INDICATIONS

- > Immobilization and support for mild spains and strains of the shoulder complex
- > Bruises and bursitis
- > Immobilization for fractures of the shoulder

FEATURES

- > Envelope design with thumb support
- > Shoulder lenght adjustable strap with foam pad for great comfort
- > Hook and loop closure for ease of adjustment
- > Ambidextrous

SIZE	ARM LENGHT (from elbow to wrist)
S	28 cm – 33 cm
M	33 cm – 38 cm
L	38 cm – 43 cm
XL	43 cm – 48 cm

INDICATIONS

- > Back pain
- > Lumbar / abdominal support

FEATURES

- > High quality breathable elastic fabric construction
- > Soft foam padding offers excellent comfort
- > 7.5 cm reinforced elastic cross-bands that can be tensioned for extra support
- > Hook and loop closure to adjust the compression to the desired intensity

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm



SOFT BACK SUPPORT

5501



BACK SUPPORT

8" 5502


**INDICATIONS**

- > Back pain
- > Lumbosacralgia
- > Muscles contractures
- > Post surgery lumbar support
- > Dorsal-lumbar disc herniation

FEATURES

- > Structure in breathable and resistant mesh elastic fabric
- > 2 removable easily shapable metal back stays integrated into the band to support the lumbar area
- > Tensionable cross elastic bands of reinforcement with hook and loop closure to adjust the compression to the desired intensity
- > Front hook and loop closure

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm

**INDICATIONS**

- > Back pain
- > Lumbar trauma
- > Muscles contractures
- > Post surgery lumbar support

FEATURES

- > Structure in breathable and resistant mesh elastic fabric
- > 4 flexible stays in the back and sides for extra lumbar support
- > Tensionable cross elastic bands of reinforcement with hook and loop closure to adjust the compression to the desired intensity
- > Front hook and loop closure

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm



SACRO LUMBAR SUPPORT

10" 5503

10"

5503



INDICATIONS

- > Back pain
- > Lumbar trauma
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Structure in breathable and resistant mesh elastic fabric
- > 4 flexible stays in the back and sides for extra support
- > Extra EVA foam pad at the back for compression and lumbar stabilization
- > Tensionable cross elastic bands of reinforcement to adjust compression and support
- > Front hook and loop closure

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm



SACRO LUMBAR SUPPORT with back pad

12"

5504



SPINAL BRACE with back pad

5505


INDICATIONS

- > Back pain
- > Lumbar trauma
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Made of breathable and resistant mesh elastic fabric for a better distribution of compression
- > 2 removable easily shapable metal back stays integrated into the band
- > EVA foam pad at the back for lumbar stabilization
- > Adjustable shoulder straps with hook and loop closure to promote good back posture
- > Tensionable cross elastic back bands to adjust lumbar compression
- > Front hook and loop closure for an easy fit

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm
XXL	112 cm – 120 cm



DELUXE SACRO LUMBAR SUPPORT 5506



INDICATIONS

- > Back pain
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Made of breathable high quality elastic fabric
- > Abdominal panels to increase abdominal compression
- > 4 flexible stays at the back for extra support
- > 3" elastic tensionable straps for more stability
- > Hook and loop closure for an easy fit

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm

INDICATIONS

- > Post surgery
- > Post pregnancy
- > Abdominal support

FEATURES

- > Made of high quality mesh elastic fabric
- > High abdominal compression
- > Highly breathable to promote the healing process
- > Front hook and loop closure to ease adjustment

SIZE	WAIST CIRCUMFERENCE
S/M	63 cm – 96 cm
L/XL	84 cm – 117 cm

ABDOMINAL SUPPORT BAND 11"

11"

5517



ABDOMINAL SUPPORT BAND 12"

12"

5508

INDICATIONS

- > Abdominal support
- > Post surgery
- > Post pregnancy

FEATURES

- > High quality durable elastic fabric
- > Made of 4 panels for a conforming fit and adjustable compression
- > Front hook and loop closure to ease adjustment
- > 12" height
- > Useful in case of particularly large surgical incisions

SIZE	WAIST CIRCUMFERENCE
S/M	63 cm – 96 cm
L/XL	84 cm – 117 cm





HYPEREXTENSION SPINE BRACE (SHORT LENGTH)

5514

HYPEREXTENSION SPINE BRACE (REGULAR LENGTH)

5515

INDICATIONS

- > Compression of spinal fractures
- > Vertebral collapse
- > Back lumbar arthrosis
- > Traumatic or pathological symptoms resulting from osteoporosis

FEATURES

- > Three-point light aluminum construction with soft foam padding for more comfort
- > Immobilization of the spine thanks to the three support points on the patient: the sternal plate on the sternum, the dorsal plate on the spine and the pelvic band on the pelvis
- > Adjustable and tilting axillary joints allow correct fitting and maximum functionality
- > Ergonomic safety closure for controlled compression

SIZE	HIPS CIRCUMFERENCE	5514 HEIGHT	5515 HEIGHT
S	60 cm – 75 cm	36 cm – 44 cm	39 cm – 47 cm
M	75 cm – 90 cm	39 cm – 47 cm	42 cm – 50 cm
L	90 cm – 105 cm	42 cm – 50 cm	45 cm – 53 cm
XL	105 cm – 115 cm	45 cm – 53 cm	48 cm – 56 cm



INDICATIONS

- > Back pain
- > Muscles contractures
- > Post surgery lumbar support
- > Abdominal support

FEATURES

- > Made of breathable high quality mesh elastic fabric
- > Abdominal panels to increase abdominal compression
- > 4 flexible stays at the back for extra support
- > 3" elastic tensionable straps for more stability
- > Hook and loop closure for an easy fit

SIZE	WAIST CIRCUMFERENCE
S	71 cm – 79 cm
M	81 cm – 89 cm
L	91 cm – 99 cm
XL	101 cm – 109 cm



BACK SUPPORT 9"

5519



INDICATIONS

- > Support for elbow instability
- > Pre and post elbow surgery

FEATURES

- > Elbow brace allows controlled range of motion for both extension and flexion
- > Detachable hand rest
- > Detachable shoulder sling
- > Ambidextrous

SIZE	FOREARM CIRCUMFERENCE
S/M	Less than 21 cm
L/XL	22 cm – 30 cm

ROM ELBOW BRACE**5309****ELBOW SUPPORT COOLMAX®****53200****INDICATIONS**

- > Supports the elbow joint area
- > Protects against bruising or scraping

FEATURES

- > Made of COOLMAX® fabric
- > Neoprene can be applied even in case of swollen or injured elbow
- > Hoop and loop closure for an easy fit
- > Ambidextrous

SIZE

ONE SIZE**INDICATIONS**

- > Immobilization and protection for wrist sprains and fractures

FEATURES

- > Anatomic plastic splints above and below the palm to provide excellent immobilization
- > The lightness of the splint together with the high quality soft foam grant an optimal comfort
- > Adjustable wrist wrap-around straps for ease of application
- > Ambidextrous

SIZE

ONE SIZE**WRIST SPLINT****53210 - left****53220 - right**

WRIST SPLINT

9"

5314 - left
5314 - right



INDICATIONS


- > Immobilization and protection of the wrist in case of sprains, tendonitis and arthrosis
- > Post fractures treatment
- > Fractures of the extremity of the radius and ulna
- > Substitute of casting

FEATURES

- > Breathable fabric
- > Adjustable hook and loop closure to offer greater stability
- > Removable aluminum splints to improve immobilization
- > Right and left

SIZE	WRIST CIRCUMFERENCE (left / right)
S	13 cm – 15 cm
M	15 cm – 17 cm
L	17 cm – 19 cm
XL	19 cm – 21 cm

BREATHABLE



INDICATIONS


- > Immobilization and protection of the wrist in case of sprains, tendonitis and arthrosis
- > Post fractures
- > Lesions of the ulnar area of the metacarpophalangeal joint of the thumb

FEATURES

- > Breathable fabric
- > Adjustable hook and loop closures to offer greater stability
- > Removable front and side aluminum splints to improve wrist and thumb immobilization
- > Right and left

SIZE	WRIST CIRCUMFERENCE (left / right)
S	13 cm – 15 cm
M	15 cm – 17 cm
L	17 cm – 19 cm
XL	19 cm – 21 cm

BREATHABLE



THUMB WITH WRIST SPLINT

5315 - left
5315 - right



THUMB WRIST SUPPORT

53030



INDICATIONS

- > Wrist and thumb immobilization and support
- > Joint instability
- > Sprains
- > Arthrosis

FEATURES

- > Made of elastic fabric with thumb support
- > Metal stay insert to keep the thumb in position
- > Adjustable wrist and thumb wrap-around straps to offer extra stability
- > Hook and loop closures for ease of application
- > Ambidextrous

SIZE	WRIST CIRCUMFERENCE
ONE SIZE	14 cm – 28 cm

INDICATIONS

- > Compression of tendon insertions in the elbow
- > Resolution of pain caused by sports injuries and by elbow overexertion (tennis elbow)
- > Prevention of epicondylitis

FEATURES

- > Anatomically shaped neoprene structure
- > Integrated pressure pad
- > Hook and loop closure for compression control
- > Elastic finishing edge

SIZE	ELBOW CIRCUMFERENCE
ONE SIZE	18 cm - 36 cm

TENNIS ELBOW SUPPORT**53070****THUMB WRIST SPLINT****53160****INDICATIONS**

- > Immobilization and support of the thumb
- > Tendinitis
- > Rhizoarthrosis
- > Post fracture

FEATURES

- > Durable and breathable perforated fabric to ensure excellent immobilization and patient comfort
- > Adjustable hook and loop closure to offer greater stability
- > Metal splint to keep thumb in place
- > Ambidextrous

SIZE	WRIST CIRCUMFERENCE
ONE SIZE	14 cm - 28 cm

INDICATIONS


- > Sprains of the thumb
- > Treatment of De Quervains syndrome

FEATURES

- > Plastic structure for excellent immobilization
- > The lightness of the brace and the soft padding grant optimal comfort
- > Adjustable hook and loop closures for easy application

SIZE	WRIST CIRCUMFERENCE
S/M	13 cm - 19 cm
L/XL	19 cm - 27 cm

THUMB WRIST SPLINT**5325-left****THUMB WRIST SPLINT****5326-right**



COOLMAX®

COOLMAX® is a breathable fabric that leaves evaporate moisture from skin making it fresh and dry and ensures maximum comfort to the patient. Coolmax® is a trademark of INVISTA.

INDICATIONS

- > Immobilization of the wrist
- > Sprains
- > Carpal tunnel syndrome
- > Tendinitis
- > Arthritis
- > Contractures

FEATURES

- > In COOLMAX® fabric
- > The neoprene structure allows it to be applied even in the event of a wrist injury
- > Adjustable hook and loop closure to facilitate adjustment

SIZE	WRIST CIRCUMFERENCE
ONE SIZE	up to 28 cm

COOLMAX® WRIST SUPPORT 53300



COOLMAX® DUAL PAD ELBOW STRAP 5345



INDICATIONS

- > Compression of tendon insertions at the elbow
- > Resolution of pain caused by sports injuries and by elbow overexertion (tennis elbow)
- > Acute and chronic epicondylitis

FEATURES

- > COOLMAX® fabric
- > Anatomically shaped neoprene
- > Silicone pads to increase pressure
- > Adjustable hook and loop closure for compression control
- > Elastic finishing edge
- > Ambidextrous

SIZE	ELBOW CIRCUMFERENCE
S	22 cm – 27 cm
M	24 cm – 29 cm
L	26 cm – 31 cm
XL	29 cm – 34 cm

INDICATIONS

- > Finger immobilization and protection
- > Fractures
- > Sprains
- > Dislocations

FEATURES

- > Plastic shell to provide superior protection and immobilization to injured finger
- > Hook and loop closures for an optimal stability
- > Lightweight and of easy application
- > Elastic finishing edge

SIZE	LENGTH OF THE SPLINT
S	5,7 cm
M	6,8 cm
L	8,3 cm
XL	9,3 cm

FINGER SPLINT 5319



SUPER ELASTIC KNEE SUPPORT

5701

INDICATIONS

- > Support and compression for weak or overworked knee
- > Patellar instability
- > Prevention in sports practice

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm



SUPER ELASTIC KNEE SUPPORT with patella opening

5702



INDICATIONS

- > Support and compression for weak or overworked knee
- > Patellar syndrome
- > Patellar tendinopathy
- > Prevention in sports practice

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > Patella opening
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm

INDICATIONS

- > Compression and patellar instability
- > Stabilization of the knee joints
- > Reduction of pain and swelling
- > Consequences of patella dislocation

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > Sides steel splints to provide additional stability
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm

SUPER ELASTIC KNEE BRACE

5703



SUPER ELASTIC KNEE BRACE with patella opening

5704



INDICATIONS

- > Compression and patellar instability
- > Stabilization of the knee joints
- > Reduction of pain and swelling
- > Consequences of patella dislocation

FEATURES

- > Durable super elastic fabric with extension in 4 ways
- > With patella opening
- > Sides steel splints to provide additional stability
- > Easy pull on application
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	30 cm – 40 cm
M	35 cm – 45 cm
L	39 cm – 49 cm
XL	43 cm – 53 cm

KNEE IMMOBILIZER 3-panels

20"

57150

INDICATIONS

- > Post-operative treatment of knee ligament reconstructions
- > Post-traumatic phase knee sprains or patella fracture
- > In emergency and first aid situations
- > To replace the casting immobilization

KNEE IMMOBILIZER 3-panels

24"

57160

FEATURES

- > Fabric with soft foam padding
- > Mouldable rigid side and back splints
- > Three adjustable panels
- > Patellar opening
- > Adjustable fastening straps with hook and loop closure
- > Ambidextrous

SIZE

ONE SIZE



ROM KNEE BRACE

19"

57110

ROM KNEE BRACE

22"

57120

INDICATIONS


- > Patella fractures
- > Post traumatic knee sprains
- > Post-operative treatment of knee ligament reconstructions
- > Early removal of casting

FEATURES

- > Aluminium hinges
- > Complete range of motion control for flexion and extension
- > Hook and loop closure straps
- > Soft polyurethane foam padding
- > Ultralight
- > Ambidextrous


SIZE	TIGHT CIRCUMFERENCE	CALF CIRCUMFERENCE
ONE SIZE	up to 70 cm	up to 48 cm





COOLPLUS®

COOLPLUS® is a revolutionary breathable fabric which, thanks to its fiber structure, leaves the skin fresh and dry.
Copyright © CHUNG SHING TEXTILE CO., LTD.



COOLPLUS® MESH HINGED KNEE BRACE

5717

INDICATIONS


- > Instability of the medial and lateral ligaments
- > Prevention and treatment of knee sprain injuries
- > Arthritic instability
- > Senile knee pathology
- > Prevention in sports activity

FEATURES

- > COOLPLUS® breathable fabric
- > Patellar opening with stabilizer
- > Front opening
- > Adjustable hook and loop closures above and below the knee to offer greater stability
- > Aluminum side splints to improve immobilization
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	33 cm – 38 cm
M	38 cm – 43 cm
L	43 cm – 48 cm
XL	48 cm – 53 cm
XXL	53 cm – 58 cm

BREATHABLE



COOLPLUS® POST-OP ROM KNEE BRACE

13"

5724

COOLPLUS® POST-OP ROM KNEE BRACE

17"

5725

INDICATIONS


- > Severe knee injuries
- > Post surgery for anterior and posterior cruciate reconstruction
- > Medial and lateral collateral ligament repair
- > Meniscal revisions

FEATURES

- > Made of COOLPLUS® breathable fabric
- > Splints with bilateral ROM mechanism
- > Patented ROM mechanism
- > Flexion-extension stop kit
- > System with extra extension stops
- > Patellar opening with stabilizer
- > Hook and loop adjustable closure straps
- > Ultralight
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
S	33 cm – 38 cm
M	38 cm – 43 cm
L	43 cm – 48 cm
XL	48 cm – 53 cm
XXL	53 cm – 58 cm

COOLPLUS® is a revolutionary reathable fabric which, thanks to its hair structure, leaves the skin fresh and dry.




ROM hinges on the lateral and medial sides to provide perfect immobilization



Patented ROM hinge allows easy and fast range of motion settings



Extra extension stop system to slow down knee extension movement and avoiding knee schock due to complete extension



Flexion / extension stop kit

Extension range settings:
0°, 15°, 30°, 45°, 60°, 75°, 90°

Flexion range settings:
0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°





COOLMAX®

COOLMAX® is a breathable fabric that leaves evaporate moisture from skin making it fresh and dry and ensures maximum comfort to the patient. Coolmax® is a trademark of INVISTA.

COOLMAX® SPORT KNEE SUPPORT 57200



INDICATIONS

- > Patellar instability
- > Consequences of patella dislocation
- > Prevention in sports practice

FEATURES

- > COOLMAX® breathable fabric
- > Adjustable hook and loop closures to offer greater stability and a comfortable fit
- > Patellar opening with stabilizer
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
ONE SIZE	up to 55 cm

INDICATIONS

- > Patellar syndrome
- > Pre-post accident patella blocking
- > Overloading of the sportsman's patellar tendon
- > Inflammations
- > Tendinopathies

FEATURES

- > COOLMAX® breathable fabric
- > Anatomically shaped orthopedic neoprene
- > With extra pad to provide patella compression
- > Adjustable hook and loop closures to offer greater stability
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
ONE SIZE	up to 55 cm

COOLMAX® KNEE SUPPORT 57300



COOLMAX® PATELLA STRAP with extra pad 57900



INDICATIONS

- > Patellar syndrome
- > Pre-post accident patella blocking
- > Overloading of the sportsman's patellar tendon
- > Inflammations
- > Tendinopathies

FEATURES

- > COOLMAX® breathable fabric
- > Anatomically shaped orthopedic neoprene
- > With extra pad to provide patella compression
- > Adjustable hook and loop closures to offer greater stability
- > Ambidextrous

SIZE	KNEE CIRCUMFERENCE
ONE SIZE	up to 60 cm

ANKLE SUPPORT

5901

INDICATIONS

- > Instability of the ankle
- > Sprains
- > Edema
- > Prevention of accidental trauma to ankle in sport
- > Useful in re-start of sporting activity

FEATURES

- > Super elastic fabric
- > Heel opening
- > Easy to wear
- > Washable
- > Elastic finishing edge
- > Ambidextrous

SIZE	ANKLE CIRCUMFERENCE
S	18 cm – 23 cm
M	21 cm – 25 cm
L	24 cm – 29 cm
XL	29 cm – 33 cm



COOLMAX® ANKLE STRAP

59200

INDICATIONS

- > From mild to moderate ankle sprains
- > Preventive bandaging

FEATURES

- > COOLMAX® breathable fabric
- > Neoprene bandaging material can be easily applied to swollen or injured ankle
- > Hook and loop closures for a comfortable and easy fit
- > One size fits all

SIZE

ONE SIZE



AIR ANKLE STIRRUP BRACE 10"

59180

AIR ANKLE STIRRUP BRACE 9"

59190

INDICATIONS

- > Prevention and treatment of ankle sprains
- > Tendinitis
- > Prevention of accidental ankle trauma in sporting activity
- > Arthrosis of the ankle joint

FEATURES

- > Structure in composite plastic material, resistant and anatomic
- > Anatomical lateral shells for ankle immobilization
- > Side air cushions to aid circulation and reduce swelling
- > Hook and loop closures for more comfort
- > Ambidextrous

SIZE

ONE SIZE

ANKLE WALKER

17"

5907

INDICATIONS

- > Acute ligament injuries
- > Rupture of the Achilles tendon
- > Malleolar fractures
- > Post-operative treatment of leg and foot fractures
- > Substitute for casting

FEATURES

- > Durable light material construction
- > Wide and soft support base to increase comfort and stability
- > EVA foot bed padding to ensure flexibility
- > The tilting sole helps promote natural stepping and reduce plantar pressure
- > Hook and loop closure straps
- > Ambidextrous

SIZE	MAX FOOT LENGTH	MAX SHOE SIZE
S	24 cm	36 - 37
M	28 cm	38 - 40
L	32 cm	41 - 44
XL	36 cm	45 - 48



ROM ANKLE WALKER

17"

5909



INDICATIONS

- > Acute ligament injuries
- > Rupture of the Achilles tendon
- > Malleolar fractures
- > Post-operative treatment of leg and foot fractures
- > Substitute for casting

FEATURES

- > Durable light material construction
- > ROM mechanism to control range of motion
- > ROM settings between 0° and 45° plantar and dorsal flexion (increments of 15°)
- > Wide and soft support base to increase comfort and stability
- > EVA foot bed padding to ensure flexibility
- > Hook and loop closure straps
- > Ambidextrous

SIZE	MAX FOOT LENGTH	MAX SHOE SIZE
S	24 cm	36 - 37
M	28 cm	38 - 40
L	32 cm	41 - 44
XL	36 cm	45 - 48

INDICATIONS

- > Foot and ankle fractures
- > Joints sprains
- > Acute knee distorsions
- > Post-operative treatment
- > Substitute for casting

FEATURES

- > Rigid construction in durable and lightweight composite material to ensure optimal support and immobilization to the foot and the ankle
- > Wide and soft support base to increase comfort and stability
- > Soft, breathable PU padding for greater comfort
- > The tilting sole helps the natural gait
- > Non-slip and wear-resistant sole for greater safety
- > Ambidextrous

SIZE	BASE LENGTH	BASE WIDTH	HEIGHT	SHOE SIZE
S	16,3 cm	7,8 cm	17,5 cm	22-26
M	20,1 cm	9,1 cm	22,5 cm	27-31
L	22,5 cm	10,4 cm	28,7 cm	32-35

PEDIATRIC WALKER

5927



ACCURACY IN ONE CLICK



DIAGNOSTICS

Diagnostics - thermometers, blood pressure monitors, pulse oximeters - **with advanced technology** and a modern, elegant and functional design, **dedicated to home and medical hospital use.**

Medi4Care THERMOMETERS are based on advanced technology that allows the whole family to measure body temperature instantaneously and accurately, at home and in hospital.


The devices have a modern, elegant and functional design.



INFRARED THERMOMETER ear / forehead

MC 31T

Ref. 410002



FEATURES

- > Multifunction: it measures the temperature in ear and forehead, displays the room temperature and acts as a clock
- > Instant measurement
- > Stores up to 10 measurements
- > Measure temperature in °C (Celsius scale) or °F (Fahrenheit scale)
- > The waterproof infrared sensor does not require the use of protection
- > Maxi backlit LCD display
- > Measurement accuracy
- > When placed in its housing, it becomes an elegant clock with the display of the room temperature
- > Auto OFF for energy saving
- > Low battery indicator
- > Eco design
- > Easy to clean



AUTO
ON/OFF



10 MEMORY
POSITIONS



MAXI BACKLIT
DISPLAY



EAR



FOREHEAD



ROOM
TEMPERATURE



TIME

Medi4Care ARM AND WRIST BLOOD PRESSURE MONITORS are indicated to monitor blood pressure both in hospital or at home. They have an innovative design with large, easy-to-read LCD display that simultaneously displays systolic and diastolic blood pressure, heart beats, date / time and the detection of IHD irregular heartbeat (Arrhythmia).

The Fuzzy Logic measurement method system ensures quick and accurate measurements according to International standards (WHO).



BLOOD PRESSURE MONITOR arm type

MC 65P

Ref. 410030



FEATURES

- > Automatic ON / OFF
- > Instant and accurate
- > Maxi multifunction display with date, time, heart rate and blood pressure
- > Indicates the room temperature on the Fahrenheit or Celsius scale
- > Low battery indicator
- > Power supply connection
- > 2 Memory positions: stores up to 60 measurements per patient and indicates the average of the last 3
- > Reminder wake-up alarm with 3 time settings
- > Pulse irregularity detector
- > Multi-information design
- > WHO indicator (World Health Organization)





BLOOD PRESSURE MONITOR arm type

MC 61P

Ref. 410020



FEATURES

- > Instant and accurate
- > Maxi multifunction display with date, time, heart rate and blood pressure
- > Low battery indicator
- > Auto ON/OFF
- > Stores up to 60 measurements and shows the average of the last 3
- > Reminder wake-up alarm with 3 time settings
- > Pulse irregularity detector
- > Multi-information design
- > WHO (World Health Organization) indicator



FEATURES

- > Instant and accurate
- > Maxi multifunction display with date, time, heart rate and blood pressure
- > Low battery indicator
- > Auto ON/OFF
- > Stores up to 90 measurements and shows the average of the last 3
- > Pulse irregularity detector
- > WHO (World Health Organization) indicator



BLOOD PRESSURE MONITOR wrist type

MC 15P

Ref. 410010



Medi4Care FINGER PULSE OXIMETERS are intended for measuring functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate to adults.

Ideal for hospital, professional, home use and for patients undergoing oxygen therapy. The color display shows the SpO₂ value, heart rate, pulse indicator, battery charge status and plethysmographic wave.



FINGERTIP PULSE OXIMETER MC 602C

Ref. 410103

FEATURES

- > Suitable for adults and children over 10 kg
- > Compact and lightweight design (60g including batteries)
- > Instant, precise, accurate and easy to use
- > LCD display adjustable in 4 different positions
- > Low battery indicator
- > Pulse rate indicator
- > Possibility of data recording up to 32 consecutive hours
- > Works with two AAA alkaline batteries (included)
- > Auto ON/OFF
- > Hanging rope for carrying around the neck



AUTO
ON/OFF



LCD
COLOR DISPLAY



FEATURES

- > Indicated for use on pediatric patients (weight less than 10 kg)
- > Compact and lightweight design (30g including battery)
- > Instant, precise, accurate and easy to use
- > Color LCD display adjustable in 4 different positions
- > Audio signal that reproduces the sound of the beat (can be deactivated)
- > Sound alarm (can be deactivated)
- > Low battery indicator
- > Pulse rate indicator
- > Possibility of data recording up to 6 consecutive hours
- > Works with a rechargeable 3,6 V button battery (power supply, USB cable and battery included)
- > Auto ON/OFF
- > With hanging rope



AUTO
ON/OFF



LCD
COLOR DISPLAY



ALARM
FUNCTION

PEDIATRIC FINGERTIP PULSE OXIMETER

MC 600

Ref. 410100

EFFICACY AND RELIABILITY IN WOUND CARE



WOUND CARE

A complete range of advanced medical devices - Polyurethane Films, Hydrocolloids, Silicone Polyurethane Foams - designed to direct clinical practice according to the most recent and proven scientific evidence, to improve the quality of care, to reduce costs and to perfect the degree of satisfaction of users and operators.

MEDI4SKIN®

SEMIPERMEABLE POLYURETHANE FILM DRESSINGS

Semi-permeable transparent PU dressings are characterized by their impermeability to water and bacteria and by their permeability to oxygen and water vapor.



Features that allow them to create and maintain the ideal moist environment for the re-epithelialization, while protecting the wound from the risk of exogenous infections.




The exudate on the wound surface, thanks to the presence of neutrophils, plays a natural antibacterial function and promotes cell growth and healing.



The perfect transparency of the film also allows proper inspection of the wound without removing the dressing.

COMFORT AND SIMPLICITY ARE EVIDENT





MEDI4SKIN®

STERILE ADHESIVE PU FILM WITH HIGH MVTR

FEATURES

- > Create the ideal moist environment for rapid healing
- > Its transparency allows continuous monitoring of the wound without removing the dressing
- > The conformability and flexibility of the film ensure maximum patient comfort
- > The rounded corners minimize the risk of accidental removal of the dressing which can thus remain in situ for longer, with obvious benefits in terms of cost-effectiveness in wound management
- > The hypoallergenic acrylic adhesive guarantees respect for the skin and prevents any irritation
- > The frame opening system makes it easy to apply the dressing
- > Equipped with an adhesive strip to note the date of application

INDICATIONS

- > As a secondary dressing
- > Superficial wounds and burns
- > Minor abrasions
- > Skin protection in the periwound area

REF	SIZE	PIECES x BOX
700107	6X7 cm	100
700112	10X12 cm	100
700120	15X20 cm	50



LOW EXUDATE




BREATHABLE



IMPERMEABLE
TO BACTERIA
AND LIQUIDS



TRANSPARENT



MEDI4SKIN[®] PLUS

STERILE ADHESIVE PU FILM ISLAND DRESSING WITH ABSORBENT NON-ADHERENT PAD

FEATURES

- > The conformability of the Polyurethane film makes the dressing adaptable to all parts of the body as a second skin
- > The impermeability of the film to liquids increases patient compliance
- > The pad retains exudate minimizing the risk of skin maceration
- > The non-adherent layer of the central pad minimizes trauma in dressing removal
- > The transparency of the film allows the operator to monitor the level of exudate present in the pad

INDICATIONS

- > Post-operative wounds
- > Lacerations
- > Abrasions

REF	SIZE	PIECES x BOX
700207	6X7 cm	50
700209	6X9 cm	50
700210	8X10 cm	50
700215	8X15 cm	50
700220	10X20 cm	25
700225	10X25 cm	25
700230	10X30 cm	25



MODERATE EXUDATE



BREATHABLE



IMPERMEABLE TO BACTERIA AND LIQUIDS



TRANSPARENT



ABSORBENT PAD



MEDI4SKIN[®] IV

STERILE ADHESIVE PU FILM DRESSING FOR THE MANAGEMENT OF INSERTION SITES

FEATURES

- > The conformability of the film and the particular shape of the dressings allow an easy, quick and aseptic application
- > The perfect transparency allows a constant monitoring of the insertion points
- > The dressing guarantees operator safety and patient comfort
- > Equipped with an adhesive strip to note the date of application and non-woven fabric strips for fixation

INDICATIONS

- > Cannulae and catheters fixation

REF	SIZE	PIECES x BOX	
700304	4,4X4,4 cm	100	
700308	6X8 cm	50	
700309	7X8,5 cm	50	
700311	8,5X10,5cm	50	



CANNULAE
FIXATION




BREATHABLE



IMPERMEABLE
TO BACTERIA
AND LIQUIDS



TRANSPARENT



MEDI4SKIN[®] ROLL

TRANSPARENT ADHESIVE PU FILM ROLL

FEATURES

- > Flexible and adaptable even to irregular areas of the body, it can be cut and used as desired
- > Ensures optimal fixation to all types of dressings and guarantees high skin breathability
- > Effectively protects against liquids and bacteria: patients can continue to bathe or shower while maintaining the dressing in situ
- > It combines the qualities of Medi4skin[®] film with great ease of use
- > The thin Polyurethane film is supported by two pre-cut films to facilitate application

INDICATIONS

- > Fixation of primary dressings
- > Minor abrasions
- > Fixation of cannulae and catheters
- > Skin protection in the friction points of drainage bags, ostomy bags, etc.

REF	SIZE	PIECES x BOX
695005	10 mX5 cm	1
695010	10 mX10 cm	1



BREATHABLE

IMPERMEABLE
TO BACTERIA
AND LIQUIDS

TRANSPARENT

MEDI4AKT[®] HD

HYDROCOLLOIDS WOUND DRESSINGS

Hydrocolloid dressings are available on the market for many years and the countless studies on them continue to confirm their effectiveness across a wide range of chronic and acute wounds.

Hydrocolloids create a moist interface between wound and dressing, facilitate autolysis and promote granulation, thus creating an environment which results optimal for wound healing.

Hydrocolloid dressings are made up of a layer of gelling material adhered to a semipermeable Polyurethane film.

The result is an absorbent and self-adhesive dressing which adheres even in humid conditions.



ALWAYS AT SERVICE OF EFFECTIVENESS

TRANSPARENT FILM


The semi-permeable Polyurethane film creates the moist environment ideal for promoting wound healing and acts as barrier to exogenous bacteria and external pollutants


The transparency of the film also allows a correct inspection of the wound without necessarily having to remove the dressing

HYDROCOLLOIDAL LAYER

The gelling layer is made of an adhesive matrix containing a combination of absorbent materials such as carboxymethylcellulose (CMC) which increase the capacity of the dressing to absorb wound exudate

The formation of a cohesive gel allows an atraumatic removal of the dressing, without damaging the newly formed tissues





MEDI4AKT® HD Thin

STERILE THIN HYDROCOLLOID DRESSING

FEATURES

- > Creates the ideal moist environment for the healing process by promoting angiogenesis, increasing the number of dermal fibroblasts and stimulating the production of granulation tissue
- > Thanks to Medi4akt® HD Thin's ability to retain moisture, the dressing helps to soften and rehydrate the necrotic tissue and slough, facilitating autolytic debridement
- > The transparency of the dressing allows to continuously monitor the wound
- > Medi4akt® HD Thin has a good absorbency and a high breathability
- > When using the dressing the absorbed exudate is transformed into a gel: this makes the removal easy and atraumatic and, by protecting nerve endings, helps reduce pain in the wound bed
- > The smooth outer coating of the dressing reduces the coefficient of friction and decreases the risk of accidental removal of the same
- > Medi4akt® HD Thin is waterproof: patients can continue to bathe or shower while keeping the dressing in place
- > Excellent cost-effectiveness ratio

INDICATIONS

- > Low to moderately exuding wounds
- > Skin ulcers
- > Pressure sores
- > Lower limb ulcers

REF	SIZE	PIECES x BOX
705010	10X10 cm	10
705015	15X15 cm	10





MEDI4AKT® HD Sacral

STERILE CONTOURED HYDROCOLLOID DRESSING

FEATURES

- > Medi4akt® HD Sacral is adhesive, waterproof and creates a barrier against microorganisms: for these reasons the contoured dressing is ideal for use in areas prone to heavy contamination, such as the sacrum
- > The area covered by Medi4akt® HD Sacral is exposed to lower pressure, stretching and friction
- > Medi4akt® HD Sacral, in the presence of devitalized tissue, facilitates autolytic debridement by creating a moist interface between wound and dressing
- > The transparency of the dressing allows to continuously monitor the wound without removing it
- > Medi4akt® HD Sacral ensures protective coverage for healthy periwound skin
- > The gel that forms when Medi4akt® HD Sacral absorbs exudate keeps the moist environment and prevents build-up of fluids on the wound
- > Easy and atraumatic removal
- > Medi4akt® HD Sacral protects nerve endings by supporting pain management
- > Excellent cost-effectiveness ratio

INDICATIONS

- > Sacral pressure sores
- > Moderately exuding wounds
- > Skin ulcers
- > Burns
- > Donor sites

REF	SIZE	PIECES x BOX
705048	15X18 cm	5



MODERATE
EXUDATE



BREATHABLE



IMPERMEABLE
TO BACTERIA
AND LIQUIDS



TRANSPARENT


MEDI4AKT[®] SF

SILICONE AND PU FOAM DRESSINGS

Polyurethane foam dressings are a tool of excellence in the management of moderately to heavily exuding wounds thanks to the high absorbency combined with the maintenance of the moist environment, high conformability, to a balanced pressure and friction protection and to the slight adherence given to the dressing by the silicone layer.

The highly breathable film acts as a barrier against infections and is water resistant.

The "cushioning" effect of the foam guarantees a balanced pressure relief.



A SOFT BREATH OF RELIEF

OUTER BREATHABLE FILM

The breathability of the outer Polyurethane film allows evaporation of excess moisture, ensuring the right moisture to the wound, avoiding maceration of the periwound area and increasing the in situ dressing time


PU FOAM

The structure of the Polyurethane foam guarantees a vertical absorption of the exudate which is retained inside the dressing thus preventing it from moving towards the surrounding tissue resulting in maceration

ADHESIVE

PERFORATED SILICONE LAYER

The uniformly perforated Silicone layer ensures a slight adherence to the periwound skin without sticking to the wound, minimizing discomfort when changing dressing





MEDI4AKT[®] SF

STERILE SILICONE AND PU FOAM DRESSINGS

FEATURES

- > It is a highly technological dressing that combines in a single product all the advantages of the absorbent pad foam and the Silicone wound contact layer
- > The high absorbency of the Polyurethane foam dressing reduces the number of changes and therefore the costs associated with wound healing
- > The absorbent cells that make up the foam are able to absorb significant quantities of exudate and to retain liquids even under compression
- > The Polyurethane film that covers the external side of the dressing prevents exudate from escaping but it is vapor permeable promoting the ideal moist environment for the healing process
- > The perforated Silicone contact layer provides optimum protection for the newly formed tissues
- > The light adherence of the Silicone layer ensures atraumatic and pain-free dressing changes
- > The conformability, flexibility and softness of the foam create a cushioning effect that reduces pressure on the wound, counteracts any frictional forces and thus guarantees comfort and protection for the patient

INDICATIONS

- > Venous ulcers
- > Diabetic ulcers
- > Skin grafts
- > Moderately to highly exuding wounds
- > Wounds with particularly fragile and delicate periwound skin

REF	SIZE	PIECES x BOX
705110	10X10 cm	5
705210	10X10 cm	10
705115	15x15 cm	5
705120	20x20 cm	5



HIGH
EXUDATE



BREATHABLE



IMPERMEABLE
TO BACTERIA
AND LIQUIDS



PERFORATED
SILICONE
LAYER



SOFT



MEDI4AKT[®] SF Sacral

STERILE CONTOURED SILICONE AND PU FOAM DRESSINGS

FEATURES

- > It is a contoured dressing that combines the characteristics of Medi4akt[®] SF with the practicality of a shape which adapts to the sacral region
- > The absorbent cells that make up the foam are able to absorb significant quantities of exudate and to retain liquids even under compression
- > The Polyurethane film that covers the external side of the dressing is breathable but waterproof and resistant to external contamination and liquids, making Medi4akt[®] SF Sacral suitable for applications in the sacral area that require a water repellent cover
- > The perforated Silicone contact layer provides optimum protection for the newly formed tissue
- > The slight adherence of the Silicone layer ensures pain-free dressing changes
- > The conformability, flexibility and softness of the foam create a cushioning effect that reduces pressure on the wound, counteracts any frictional forces and thus guarantees comfort and protection for the patient

INDICATIONS

- > Pressure ulcers in the sacral area
- > Diabetic ulcers
- > Skin grafts
- > Moderately to highly exuding wounds
- > Wounds with particularly fragile and delicate periwound skin

REF	SIZE	PIECES x BOX
705128	15X18 cm	5



HIGH
EXUDATE



BREATHABLE



IMPERMEABLE
TO BACTERIA
AND LIQUIDS



PERFORATED
SILICONE
LAYER



SOFT



MEDI4AKT®*F Tracheo*

STERILE TRACHEOSTOMY AND DRAINAGE PU FOAM DRESSINGS

FEATURES

- > The absorbed liquids are evenly distributed inside the foam which conforms perfectly to the tracheostomy or drainage tube ensuring maximum patient comfort
- > The high absorbency of the Polyurethane foam dressing reduces the number of changes and the costs associated with wound management
- > Key-hole ∅ 12-20mm

INDICATIONS

- > Management of exudates around the graft sites of the tracheostomy or drainage tube

REF	SIZE	PIECES x BOX
705310	10X10 cm	10



MODERATE EXUDATE



BREATHABLE



KEY-HOLE



IMPERMEABLE TO BACTERIA AND LIQUIDS



SOFT

A COMPLETE RANGE OF
ADVANCED MEDICAL DEVICES
DESIGNED TO DIRECT CLINICAL
PRACTICE ACCORDING TO THE MOST
RECENT AND PROVEN SCIENTIFIC
EVIDENCE, TO IMPROVE THE QUALITY
OF CARE AND TO PERFECT THE
DEGREE OF SATISFACTION OF
USERS AND OPERATORS.





UNIONE EUROPEA



REPUBBLICA ITALIANA



REGIONE SICILIANA

P  **FESR**
SICILIA 2014-2020